



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

ONYX

MODEL: [RE-ONYX]



DESCRIPTION:

Deep stylistic curves and lenses combine style with lateral x-ray radiation protection. The Onyx frame features a ventilated nose bridge providing grip, reduced fogging, and comfort.

DESIGN:

Frame style: Wraparounds
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 71.8g
Frame size: 68-19-135
Side Shields: None

PRESCRIPTION:

Options: N/A
Range: N/A

Minimum Progressive PD: N/A
Base Curve: 8

FEATURES:

High curve lenses, ventilated nose bridge

COLORS:

Shiny White (SWHT), Matte Black (MBLK)



SWHT



MBLK