



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

PANORAMIC FACE SHIELD



MODEL: [400PSM]

DESCRIPTION:

The Panoramic Face Shield utilizes a lead impregnated acrylic shield suspended from a dual adjustable headpiece. This lead face shield provides additional x-ray protection around the sides of the head while sacrificing some front protection. Because of the multi-adjustable headpiece, you can angle the mask to protect from radiation in many different positions as well as adjust the size to comfortably fit your head. *Mask should be worn with leaded eyewear protection.

DESIGN:

Model: 400PSM

Frame style: Face Shield

Front lens: 0.10mm Lead Acrylic

Frame Material: Plastic

Weight: 350g

Frame size: Adjustable Head Piece

PRESCRIPTION:

Options: N/A

Range: N/A

Minimum Progressive PD: N/A

Base Curve: N/A

