



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# PETITE METAL

MODEL: [ RE-007 ]

## DESCRIPTION:

The model 007 is a classic style made of optical steel with bronze electroplate. This durable frame comes standard with 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

## DESIGN:

Frame style: Metals  
Front lens: 0.75mm Pb Glass  
Frame Material: Optical Steel

Weight: 63g  
Frame size: 52-22-150  
Side Shields: 0.50mm (included)

## PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive  
Range: +8.00 to -9.00

Minimum Progressive PD: 58  
Base Curve: 4-6

## FEATURES:

Adjustable silicone nose pads, 0.50mm side shields included, versatile fit for all face types, lightweight metal frame

## COLORS:

Gold (GLD), Light Gunmetal (GUNM)



GUNM



GLD

