



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# RAY BAN 2016 DADDY-O



MODEL: [ RE-RB-2016 ]

**DESCRIPTION:** Ray Ban's Daddy-O lead glasses are a sporty frame that boast optimal coverage for maximum x-ray protection.

**DESIGN:** Frame style: Brand name - Wraparound fit  
Front lens: 0.75mm Pb Glass  
Frame Material: Nylon  
Weight: 61g  
Frame size: 59-17-125  
Side Shields: 0.50mm Pb Strips (optional)

**PRESCRIPTION:** Options: Single Vision, Bifocal, Progressive  
Range: +8.00 to -9.00  
Minimum Progressive PD: 65  
Base Curve: 6

**FEATURES:** Lightweight frame, excellent wraparound coverage for optimal x-ray protection

**COLORS:** Matte Black (MBLK)

