



**PROTECH**  
Medical

**PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# RAY BAN 3534

MODEL: [ RE-RB-3534 ]

## DESCRIPTION:

Ray Ban's 3534 lead glasses boast a sporty, masculine frame perfect for everyday use. Brushed aluminum temples, adjustable nose bridge, and rubberized ear socks ensure your eyewear are perfectly balanced and have a personalized fit.

## DESIGN:

**Frame style:** Brand name - Metals  
**Front lens:** 0.75mm Pb Glass  
**Frame Material:** Metal

**Weight:** 63g  
**Frame size:** 59-20-135  
**Side Shields:** Not available

## PRESCRIPTION:

**Options:** None  
**Range:** Not applicable

**Minimum Progressive PD:** Not applicable  
**Base Curve:** 6

## FEATURES:

Rubberized ear socks for added comfort and grip, adjustable silicone nose pads for personalized fit

## COLORS:

Black (BLK), Gunmetal (GUNM)



BLK

*Ray-Ban*



GUNM