



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

RAY BAN 4015

MODEL: [RE-RB-4105]

DESCRIPTION:

Ray Ban's 4105 wayfarer lead glasses have been constructed into a compact, folding style to be easily portable while still being highly durable.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Nylon

Weight: 62g
Frame size: 50-23-140
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 65
Base Curve: 6

FEATURES:

Lightweight frame, excellent wraparound coverage for optimal x-ray protection

COLORS:

Matte Black (MBLK)



Ray-Ban

