



PROTECH
Medical

PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone:	561.627.9769
Email(USA):	sales@protechmed.com
Email(Intl.):	global@protechmed.com
Email(LatAm.):	brazil@protechmed.com

RAY BAN 4165F JUSTIN

MODEL: [RE-RB-4165]

DESCRIPTION:

Inspired by the iconic Wayfarer sunglasses, Ray Ban's 4165 Justin lead glasses are redesigned with slightly larger rectangle lenses and a rubberized frame finish to really give it that extra edge.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Nylon

Weight: 73g
Frame size: 55-18-145
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 65
Base Curve: 6

FEATURES:

Rubberized frame finish for no-slip grip and comfort, durable nylon frame

COLORS:

Matte Black (MBLK)



Ray-Ban

