



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

RAY BAN 4181

MODEL: [RE-RB-4181]

DESCRIPTION: Designed to fit all face types, Ray Ban's 4181 lead glasses feature a larger plastic frame that is both comfortable and stylish.

DESIGN: **Frame style:** Brand name - Classic fit **Weight:** 71.4g
Front lens: 0.75mm Pb Glass **Frame size:** 57-16-145
Frame Material: Injected plastic **Side Shields:** Not available

PRESCRIPTION: **Options:** Single Vision, Bifocal, Progressive **Minimum Progressive PD:** 65
Range: +8.00 to -9.00 **Base Curve:** 6

FEATURES: Versatile fit for all face types (unisex)

COLORS: Tortoise (TORT)

