



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

RAZER

MODEL: [RE-RAZER]



DESCRIPTION:

The Razer's blends a sleek design into a lightweight and durable frame. Its premium features such as a ventilated nose bridge, rubberized temples, and wide-angle lenses maximize comfort and protection from x-ray. The modern style and use of light-weight materials make the Razer an excellent choice for long-procedures.

DESIGN:

Frame style: Wraparounds
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 62.5g
Frame size: 58-17-135
Side Shields: 0.50mm Pb strips (optional)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00

Minimum Progressive PD: 59
Base Curve: 6

FEATURES:

Ventilated Nose Bridge, Rubberized Temples for Grip and Comfort, Modern Matte Color Finishes, High Rx Options

COLORS:

Matte Green (**MLGN**), Matte Black (**MBLK**), Matte Gray (**MGRY**), Matte Blue (**MBLU**), Matte Red (**MRED**)



MLGN



MBLK



MBLU



MRED



MGRY



Clear Side Shield Insert