



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

TUGGA BABY

MODEL: [RE-BABY]

DESCRIPTION:

Protech's Tugga Baby pediatric glasses provide essential protection during medical fluoroscopy and feature a universal fit with its adjustable headband.

DESIGN:

Frame style: Goggles

Front lens: 0.75mm Pb Glass

Frame Material: Polycarbonate

Weight: 42g

Frame size: 46-18-adj

Side Shields: None

PRESCRIPTION:

Options: None

Range: Not Applicable

Minimum Progressive PD: Not Applicable

Base Curve: 4

FEATURES:

Adjustable head band for universal fit

COLORS:

Clear Pink (PNK), Clear Blue (BLU), Clear Green (GRN)



PNK



BLU



GRN