



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# WILEY-X AIRRAGE

MODEL: [ RE-WX-AIRR ]

## DESCRIPTION:

Designed to withstand high impact, the Wiley-X Airrage lead glasses feature a facial foam cavity seal and high curved lenses for maximum comfort and protection against x-ray

## DESIGN:

**Frame style:** Brand name - Wraparound fit  
**Front lens:** 0.75mm Pb Glass  
**Frame Material:** Polycarbonate

**Weight:** 67.5g  
**Frame size:** 61-18-124  
**Side Shields:** Not available

## PRESCRIPTION:

**Options:** None  
**Range:** Not applicable

**Minimum Progressive PD:** Not applicable  
**Base Curve:** 8

## FEATURES:

Highly durable frame, facial foam cavity for added comfort, highly curved lenses for optimal x-ray protection

## COLORS:

Matte Black (MBLK)

