



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

WILEY-X BOSS

MODEL: [RE-WX-BOSS]

DESCRIPTION:

The Wiley-X Boss lead glasses rest comfortably close to your face while the removable Facial Cavity Seal prevents all small dust and dirt particles from getting into your eyes. These protective eyewear boast highly curved lenses for excellent x-ray protection.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 84.5g
Frame size: 68-18-125
Side Shields: Not available

PRESCRIPTION:

Options: None
Range: Not applicable

Minimum Progressive PD: Not applicable
Base Curve: 8

FEATURES:

Highly durable frame, facial foam cavity for added comfort, highly curved lenses for optimal x-ray protection, maximum splash protection

COLORS:

Matte Black (MBLK)

