



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

WILEY-X BRICK

MODEL: [RE-WX-BRICK]

DESCRIPTION:

The Wiley-X Brick lead glasses offer a flexible rubber temple for support during intense activities while the removable Facial Cavity Seal provides maximum splash protection. These protective eyewear boast highly curved lenses for excellent x-ray protection.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 74g
Frame size: 63-18-120
Side Shields: Not available

PRESCRIPTION:

Options: None
Range: Not applicable

Minimum Progressive PD: Not applicable
Base Curve: 8

FEATURES:

Highly durable frame, rubberized flexible temple arms, removable facial foam cavity for added comfort, highly curved lenses for optimal x-ray protection

COLORS:

Matte Black (MBLK)

