



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

WILEY-X CONTOUR

MODEL: [RE-WX-CON]

DESCRIPTION:

Lightweight and comfortable, the Wiley-X Contour lead glasses feature removable side shield inserts to ensure full lateral protection from x-ray.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Nylon

Weight: 62.3g
Frame size: 54-17-140
Side Shields: 0.50mm Pb Inserts (included)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4

FEATURES:

Highly durable frame, removable 0.50mm side shield inserts

COLORS:

Gloss Black and Brown Stripe (BRWN), Matte Black (MBLK)



BRWN



MBLK

