



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

WILEY-X P17

MODEL: [RE-WX-P17]



DESCRIPTION:

The Wiley-X P-17 lead glasses are low-profile and lightweight, making them a fan favorite. These protective frames boast highly curved lenses for maximum x-ray protection. The P-17 eyewear are ideal for medium to large facial structures.

DESIGN:

Frame style: Brand name - Wraparound Fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 68g
Frame size: 60-18-125
Side Shields: Not available

PRESCRIPTION:

Options: None
Range: Not applicable

Minimum Progressive PD: Not applicable
Base Curve: 6-8

FEATURES:

Highly durable frame, highly curved lenses for optimal x-ray protection, ideal for medium to large face types

COLORS:

Shiny Black (SBLK), Matte Black (MBLK), Brown (BRWN)



MBLK



SBLK



BRWN