



**PROTECH**  
Medical

**PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# XGUARD

MODEL: [ RE-XGUARD ]



## DESCRIPTION:

Tired of the basic goggle selections for radiation protection? So are we! The X Guard goggle features a premium, clear insert, lead strip side shields for lateral coverage, and tasteful color options.

## DESIGN:

**Frame style:** Goggles

**Front lens:** 0.75mm Pb Glass

**Frame Material:** Polycarbonate

**Weight:** 81.8g

**Frame size:** 58-18-135

**Side Shields:** 0.50mm Pb strips (optional)

## PRESCRIPTION:

**Options:** Single Vision, Bifocal, Progressive

**Range:** +8.00 to -10.00

**Minimum Progressive PD:** 62

**Base Curve:** 4

## FEATURES:

Wide field of view, removable side shield insert

## COLORS:

Black/Woodgrain (**WOOD**), Matte Black (**MBLK**)



WOOD



Removable side shield insert



MBLK