

APRON SIZE CHART (MALE)

GENERAL APRON SIZING


	XS	Small	Med	Large	XL	2X	3X
Height	5'2"	5'4"	5'8"	5'11"	6'2"+		
Chest	39"	42"	45"	48"	53"	57"	60"
Waist	32	34	38	42	46	50	54
Hips	34	36	40	44	48	52	56
Front Length	34	34	37	40.5	42.5	42.5	45.5
Front Width	19	21	22	24	26	28	30
Vest Length	22	22	24	25	25	25	28
Skirt Length	20	20	20	21	21	22	22
WABR Length	36	36	38	40.5	42.5	42.5	45.5

TRI-TAB APRON SIZING

	Small	Med	Large	XL
Chest	39"	42"	46"	49"
Length	25"	27"	27.5"	28.5"

APRON SLEEVE SIZING

	Small	Med	Large	XL	2XL
Length (Regular)	11"	11"	11"	11"	11"
Length (Short)	8.5"	8.5"	8.5"	8.5"	8.5"
Circumference	15"	16-17"	18-19"	19-20"	19-20"



KIDS*

	Small	Medium	Large
Dimensions	12 x 24"	14 x 26"	16 x 28"
Ages	Ages 1-3	Ages 4-6	Ages 7-10

*Velcro Adjustable & Tie Aprons

Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Looking for a made-to-measure apron? Contact us and provide your measurements for a fully custom apron!



APRON SIZE CHART (FEMALE)

GENERAL APRON SIZING

	XS	Small	Med	Large	XL	2X	3X
Height	5'1"	5'2"	5'5"	5'8"	5'10"	5'11"	
Chest	36"	38"	40"	43"	46"	49"	52"
Waist	25	27	29	32	35	40	44
Hips	35	37	39	42	46	48	52
Front Length	33	34	36	38.5	40.5	42.5	42.5
Front Width	18	19	20	21	23	24	26
Vest Length	19	22	22	22	22	23	23
Skirt Length	20	20	22	22.5	22.5	22.5	24
WABR Length	33	34	37	40.5	42.5	42.5	44.5

TRI-TAB APRON SIZING

	Small	Med	Large	XL
Chest	35"	37"	41"	44"
Length	22"	22"	22.5"	23"


APRON SLEEVE SIZING

	Small	Med	Large	XL	2XL
Length (Regular)	11"	11"	11"	11"	11"
Length (Short)	8.5"	8.5"	8.5"	8.5"	8.5"
Circumference	15"	15"	16-17"	16-17"	17-18"

KIDS*

	Small	Medium	Large
Dimensions	12 x 24"	14 x 26"	16 x 28"
Ages	Ages 1-3	Ages 4-6	Ages 7-10

*Velcro Adjustable & Tie Aprons



Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Looking for a made-to-measure apron? Contact us and provide your measurements for a fully custom apron!



APRON MEASUREMENTS

CHEST MEASUREMENT



Measure around the largest circumference of the bust area, keeping the tape parallel to the floor with the client's arms at their sides.

WAIST MEASUREMENT




Measure around the smallest circumference of the torso. If you measure at waist level it may be inaccurate.

HIP MEASUREMENT




Measure around the largest circumference of the hips and buttocks with the tape parallel to the floor: generally 8-9 inches down from the waist.

FRONT LENGTH MEASUREMENT



Measure from the top of the shoulder over the bust to the desired length. Do not contour the tape measure in under the bust.

VEST LENGTH MEASUREMENT



Measure from the top of the shoulder over the bust to the desired length. Do not contour the tape measure under the bust.

SKIRT LENGTH MEASUREMENT



Measure from the top of the pant line (waist) and over the outside of the hip down to the desired length.



APEX ADJUSTMENT

Measure from the top of the shoulder (towards the neck) over the bust to the nipple. The apex measurement determines whether adjustments need to be made to the shoulder area. This helps to avoid shoulder gapping or tight arm holes in your apron. For example, a male large height of 5'10" may be fitted into a 3X due to their measurements. This would likely result in excessive shoulder material. Therefore, an apex adjustment of -1" would be applied.