

APRON SIZE CHART (FEMALE)

- Sizing is in inches -

GENERAL APRON SIZING

	XS	Small	Med	Large	XL	2XL	3XL
Height	5'1"	5'2"	5'5"	5'8"	5'10"	5'11"	
Chest	36	38	40	43	46	49	52
Waist	25	27	29	32	35	40	44
Hips	35	37	39	42	46	48	52
Frontal Length	33	34	36	38.5	40.5	42.5	42.5
Frontal Width	18	19	20	21	23	24	26
Vest Length	19	22	22	22	22	23	23
Skirt Length	20	20	22	22.5	22.5	22.5	24
WABR Length	33	34	37	40.5	42.5	42.5	44.5

TRI-TAB VEST SIZING

	Small	Med	Large	XL	2XL
Chest	38	40	43	46	49
Vest Length	22"	22"	22.5"	23"	26"

APRON SLEEVE SIZING

	Small	Med	Large	XL	2XL
Length (Regular)	11"	11"	11"	11"	11"
Length (Short)	8.5"	8.5"	8.5"	8.5"	8.5"
Circumference	15"	15"	16-17"	16-17"	17-18"

KIDS*

	Small	Medium	Large
Dimensions	12 x 24"	14 x 26"	16 x 28"
Ages	Ages 1-3	Ages 4-6	Ages 7-10

*Velcro Adjustable & Tie Aprons

General Instructions

Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Standard sizes shown above apply as upper limits. If your measurements fall above a size's specified limit, we recommend you size up for additional allowance, aka "wiggle room." There is NO standard size allowance added so we recommend that you incorporate an additional 1-2" of "wiggle room" for the Chest measurement.

Custom Sizing Allowances aka "Wiggle Room"


- All Styles XS to XL, by default, will be made with a +1" allowance in the Chest + Skirt waist - aka "wiggle room."
- All Styles 2XL to 4XL, by default, will be made with a +2" allowance in the Chest + Skirt waist - aka "wiggle room."
- If you do not want "wiggle room", you must specify.

Quality Assurance Size Tolerances

For all apron styles, our allowed size tolerances are -0 (no underage) up to +1".

Length Adjustments

Length adjustments vary depending on your height and the required length adjustment to ensure your apron reaches around your knee cap. Any additional inches (up to 6") incur extra fees.



APRON SIZE CHART (FEMALE)

- Sizing is in centimeters -

GENERAL APRON SIZING

	XS	Small	Med	Large	XL	2XL	3XL
Height	154.94	157.48	165.1	172.72	177.8	180.34	
Chest	91.4	96.5	101.6	109.2	116.8	124.5	132
Waist	63.5	68.6	73.7	81.3	88.9	101.6	111.8
Hips	88.9	94	99.1	106.7	116.8	121.9	132.1
Frontal Length	83.8	86.4	91.4	97.8	102.9	108	108
Frontal Width	45.7	48.3	50.8	53.3	58.4	61	66
Vest Length	48.3	55.9	55.9	55.9	55.9	58.4	58.4
Skirt Length	50.8	50.8	55.9	57.1	57.1	57.2	61
WABR Length	83.8	86.4	94	102.9	108	108	113

TRI-TAB VEST SIZING

	Small	Med	Large	XL	2XL
Chest	96.5	101.6	109.2	116.8	124.5
Vest Length	55.9	55.9	57.1	58.4	66

APRON SLEEVE SIZING

	Small	Med	Large	XL	2XL
Length (Regular)	27.9	27.9	27.9	27.9	27.9
Length (Short)	21.6	21.6	21.6	21.6	21.6
Circumference	38.1	38.1	40.6-43.2	40.6-43.2	43.2-45.7

KIDS*

	Small	Medium	Large
Dimensions	30.5 x 60.9	35.5 x 66	40.6 x 71.1
Ages	Ages 1-3	Ages 4-6	Ages 7-10

*Velcro Adjustable & Tie Aprons

General Instructions

Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Standard sizes shown above apply as upper limits. If your measurements fall above a size's specified limit, we recommend you size up for additional allowance, aka "wiggle room." There is NO standard size allowance added so we recommend that you incorporate an additional 2.5-5cm of "wiggle room" for the Chest measurement.

Custom Sizing Allowances aka "Wiggle Room"

- All Styles XS to XL, by default, will be made with a +2.5cm allowance in the Chest + Skirt waist - aka "wiggle room."
- All Styles 2XL to 4XL, by default, will be made with a +5cm allowance in the Chest + Skirt waist - aka "wiggle room."
- If you do not want "wiggle room", you must specify.

Quality Assurance Size Tolerances

For all apron styles, our allowed size tolerances are -0 (no underage) up to +2.5cm.

Length Adjustments

Length adjustments vary depending on your height and the required length adjustment to ensure your apron reaches around your knee cap. Any additional centimeters (up to 15.2cm) incur extra fees.

