

FULL PROTECTION

APRON FEATURES

- +** **Adjustable Buckle Set**
Adjustable buckle closures ensure a secure fit.
- +** **Back Relief Belt**
Elastic lumbar support belt with buckle closure system.
- +** **Seam-to-seam Protection**
Seam-to-seam, standardized protection in Front (no change in protection at overlap)

- Shoulder Padding**
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro and Buckle Closures**

ADDITIONAL CUSTOMIZATION OPTIONS

- | | | |
|----------------|-------|--------------|
| Thyroid Collar | Wings | Extra Pocket |
| Embroidery | Logo | Sleeve |

Core Material Weight (0.25mm) Thickness Options

Pb	Lead	3400 gsm	0.25, 0.50mm
☐	Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☐	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

BALANCE™ WRAPAROUND BACK RELIEF APRON

MODEL: [AP-WABR]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wraparound Back Relief apron features velcro and buckle closures to ensure a proper fit. The Balance configuration features **seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration).



Example of a 0.50mm LE Balance Wraparound Apron Lead Configuration

