

FULL PROTECTION

APRON FEATURES

+ Adjustable Buckle Set
Adjustable buckle closures ensure a secure fit.

+ Back Relief Belt
Elastic lumbar support belt with buckle closure system.

+ Seam-to-seam Protection
Seam-to-seam, standardized protection in Front (no change in protection at overlap)

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro and Buckle Closures

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket
Embroidery	Logo	Sleeve

Core Material	Weight (0.25mm)	Thickness Options
---------------	-----------------	-------------------

 Lead	3500 gsm	0.25, 0.50mm
 Prolite	2850 gsm	0.25, 0.35, 0.50mm
 Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
 Premier	2600 gsm	0.25, 0.35, 0.50mm
 Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

BALANCE™ WRAPAROUND BACK RELIEF APRON

MODEL: [AP-WABR]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wraparound Back Relief apron features velcro and buckle closures to ensure a proper fit. The Balance configuration **features seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration).



Example of a 0.50mm LE Balance Wraparound Apron Lead Configuration

