

## FRONTAL PROTECTION

# APRON FEATURES

**+** Anti-slip fabric (interior)  
Reduces slippage and improves fit.

**+** Built-in Thyroid Collar  
A built-in, velcro thyroid collar provides added protection for your thyroid

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**No Closure**

Available in adult and pediatric sizes

## ADDITIONAL CUSTOMIZATION OPTIONS

Embroidery      Logo

**Core Material**      **Weight (0.25mm)**      **Thickness Options**

<b>Pb</b>	Lead	3400 gsm	0.25, 0.50mm
<b>PL</b>	Prolite	2850 gsm	0.25, 0.35, 0.50mm
<b>PR</b>	Premier	2600 gsm	0.25, 0.35, 0.50mm
<b>N</b>	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# DENTAL APRON

MODEL: [ AP-DA ]

Features a non-slip microfiber material that ensures the apron lays in place over the patient. There are no closures on the dental apron. A built-in thyroid collar is included.

**Adult sizing: 28"W x 31"L; Pediatric sizing: 24"W x 23"L**



## FRONTAL PROTECTION

# APRON FEATURES

**+** Anti-slip fabric (interior)  
Reduces slippage and improves fit.

**+** Built-in Thyroid Collar  
A built-in, velcro thyroid collar provides added protection for your thyroid

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**No Closure**

Available in adult and pediatric sizes

## ADDITIONAL CUSTOMIZATION OPTIONS

Embroidery      Logo

Core Material      Weight (0.25mm)      Thickness Options

Pb	Lead	3400 gsm	0.25, 0.50mm
☘	Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☘	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# PEDIATRIC DENTAL APRON

MODEL: [ AP-DA-CH ]

Features a non-slip microfiber material that ensures the apron lays in place over the patient. There are no closures on the dental apron. A built-in thyroid collar is included.

**Adult sizing: 28"W x 31"L; Pediatric sizing: 24"W x 23"L**

