

EZ VEST APRON

MODEL: [AP-VEZ]

A more efficient and ergonomic vest design providing EZ comfort! The EZ Vest features a single front panel and dual back panels with large, criss-crossing, elastic strap closure tabs. The shortened back panels enable better closure maneuverability and wearability while the larger elastic straps promote a secure fit and reduce back and shoulder fatigue.



APRON FEATURES

Full Protection Apron

IEC 60% Frontal Design

Compliant with IEC 60% frontal design requirement

Elastic Straps for Support

Provides enhanced adjustability and lumbar support.

Single Front Panel

Eliminates excess overlapping panels achieving a 10% weight reduction (compared to other wraparound and VS models).

Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

Additional Customization Options

- Thyroid Collar
- Embroidery
- Logo
- Extra Pocket
- Apron Sleeve
- Axillary Wings
- Drop Shoulders

	Core Material	Weight (0.25mm)	Thickness Options
	Lead	3500 gsm	0.25, 0.50mm
	Prolite	2850 gsm	0.25, 0.35, 0.50mm
	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
	Premier	2600 gsm	0.25, 0.35, 0.50mm
	Nano	2300 gsm	0.35, 0.50mm



For color options, see pages refer to apron fabric color page.



*CE applies to Lead and Prolite IEC from 0.25mm to 1.00mm.

EZ SKIRT APRON

MODEL: [AP-SKEZ]

The EZ Skirt features a built-in elastic belt which provides EZ security and comfort! Made with a single back panel and 2 wider front panels, the EZ skirt's ergonomic design maximizes wearability and protection. Anti-slip fabric is included on the built-in elastic belt.



APRON FEATURES

Full Protection Apron

IEC 60% Frontal Design

Compliant with IEC 60% frontal design requirement

Integrated Elastic Belt

Provides enhanced adjustability and lumbar support. EZ belt system makes skirt easy to don.

Contoured Back Panel

Curved skirt panel allows for easier movement.

Non-slip Backing






Reduces slippage and improves fit.

Hanging Tabs

Velcro and Belt Closure

Additional Customization Options

- Embroidery
- Logo
- Skirt Pocket
- Suspenders

Core Material		Weight (0.25mm)	Thickness Options
	Lead	3500 gsm	0.25, 0.50mm
	Prolite	2850 gsm	0.25, 0.35, 0.50mm
	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
	Premier	2600 gsm	0.25, 0.35, 0.50mm
	Nano	2300 gsm	0.35, 0.50mm



For color options, see pages refer to apron fabric color page.



*CE applies to Lead and Prolite IEC from 0.25mm to 1.00mm.