

FULL PROTECTION

## APRON FEATURES

**+** Adjustable Buckle Set  
Adjustable buckle closures ensure a secure fit.

**+** Back Relief Belt  
Elastic lumbar support belt with buckle closure system.

**+** Extended Front Panels  
Front panels are extended to 3/4 overlap coverage where LE protection is doubled

Shoulder Padding  
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro and Buckle Closures



# ESSENTIAL™ WRAPAROUND BACK RELIEF APRON

MODEL: [ AP-WAQ ]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wraparound Back Relief apron features velcro and buckle closures to ensure a proper fit. The Essential apron configuration **features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection. The Essential Wraparound is available in standard sizing only.** For a made-to-measure Wraparound apron, see Balance Wraparound.

### ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket
Embroidery	Logo	Sleeve

Core Material	Weight (0.25mm)	Thickness Options
---------------	-----------------	-------------------

Core Material	Weight (0.25mm)	Thickness Options
Lead	3500 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
Premier	2600 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.



**Example of a 0.50mm LE Essential Wraparound Apron Lead Configuration**

