

FULL PROTECTION

APRON FEATURES

- +** **Elastic Back**
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- +** **Single Front Panel**
Eliminates excess overlapping panels achieving a 10% weight reduction (compared to VS model).
- +** **Tri-tab Closure Design**
Ensures a proper fit.
- +** **Non-slip Backing (skirt)**
Reduces slippage and improves fit.
- +** **Shoulder Padding**
Premium foam padding in the shoulders provides extra comfort.
- +** **Apron Pocket**
- +** **Multi-tab Velcro Closure**
- +** **Hanging Tabs on Skirt**

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar Extra Pocket Embroidery
Logo Sleeve Suspenders (skirt)



MODEL: [AP-TVEST] (vest only) [AP-SKIRT] (skirt only) [AP-TVS] (vest and skirt combo)

TRI-TAB VEST AND SKIRT APRON

With an easy-on, easy-off multi tab (velcro) closure design, the Tri-tab is a popular style for vest and skirt combinations. Velcro closures are located on the left shoulder and around the hips. An elastic back allows for a flexible fit while providing premium comfort and back support. The Tri-tab Vest Skirt reduces excess overlapping lead and is **up to 10% lighter than traditional vest models.**



Core Material Weight (0.25mm) Thickness Options

Pb	Lead	3500 gsm	0.25, 0.50mm
P	Prolite	2850 gsm	0.25, 0.35, 0.50mm
IEC	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
N	Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.