

FULL PROTECTION

## APRON FEATURES

- +** **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- +** **Single Front Panel**  
Eliminates excess overlapping panels achieving a 10% weight reduction (compared to VS model).
- +** **Tri-tab Closure Design**  
Ensures a proper fit.
- +** **Non-slip Backing (skirt)**  
Reduces slippage and improves fit.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Multi-tab Velcro Closure**

**Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Extra Pocket	Embroidery
Logo	Sleeve	Suspenders (skirt)

Core Material	Weight (0.25mm)	Thickness Options
---------------	-----------------	-------------------

Core Material	Weight (0.25mm)	Thickness Options
Pb Lead	3400 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.25, 0.35, 0.50mm



# TRI-TAB VEST AND SKIRT APRON

MODEL: [AP-TVEST] (vest only) [AP-SKIRT] (skirt only) [AP-TVS] (vest and skirt combo)

With an easy-on, easy-off multi tab (velcro) closure design, the Tri-tab is a popular style for vest and skirt combinations. Velcro closures are located on the left shoulder and around the hips. An elastic back allows for a flexible fit while providing premium comfort and back support. The Tri-tab Vest Skirt reduces excess overlapping lead and is **up to 10% lighter than traditional vest models.**



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.