

## FRONT PROTECTION

# APRON FEATURES



### Wide Body Back Panels

The comfortable panels secure your apron and provides lumbar support.

### Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

### Apron Pocket

### Velcro Closure

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Wings
Extra Pocket	Embroidery	Backstrap
Logo	Sleeve	Back Mesh

## Core Material Weight (0.25mm) Thickness Options

Pb	Lead	3500 gsm	0.25, 0.50mm
P	Prolite	2850 gsm	0.25, 0.35, 0.50mm
IEC	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
N	Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

# VELCRO ADJUSTABLE APRON

MODEL: [ AP-VA ]

Do you need comfort and adjustability at an affordable price? The Velcro Adjustable apron has crisscrossing back straps that provide back support while securing the apron around your waist. Wide back straps prevent the apron from slipping around your shoulder and back area. Pediatric sizing available.

