

FRONT PROTECTION

APRON FEATURES



Wide Body Back Panels
The comfortable panels secure your apron and provides lumbar support.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Wings
Extra Pocket	Embroidery	Backstrap
Logo	Sleeve	Back Mesh

Core Material Weight (0.25mm) Thickness Options

Core Material	Weight (0.25mm)	Thickness Options
Pb Lead	3400 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

VELCRO ADJUSTABLE APRON

MODEL: [AP-VA]

Do you need comfort and adjustability at an affordable price? The Velcro Adjustable apron has crisscrossing back straps that provide back support while securing the apron around your waist. Wide back straps prevent the apron from slipping around your shoulder and back area. Pediatric sizing available.

