

FULL PROTECTION

APRON FEATURES



Elastic Back

Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.



Single Front Panel

Eliminates excess overlapping panels achieving a 15% weight reduction (compared to other wraparound and VS models).



Tri-tab Closure

Ensures a proper fit.

Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar Wings Extra Pocket
Embroidery Logo Sleeve

Core Material Weight (0.25mm) Thickness Options

Pb	Lead	3500 gsm	0.25, 0.50mm
☐	Prolite	2850 gsm	0.25, 0.35, 0.50mm
☐	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☐	Nano	2300 gsm	0.35, 0.50mm



WRAPAROUND FLEXBACK APRON

MODEL: [AP-WAFB]

Want the ultimate in comfort and protection? The Wraparound Flexback apron combines the best qualities of our flexback apron with the full protection and back support of our wraparound back relief apron. The elastic stretch-back shifts weight from your shoulders to your hips while the tri-tab velcro closure system provides convenient, easy-on easy-off wearability. The Wrap Around Flexback eliminates excess overlapping lead and is up to **15% lighter than traditional wrap styles**.



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.