

FULL PROTECTION

APRON FEATURES

- +** **Elastic Back**
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- +** **Single Front Panel**
Eliminates excess overlapping panels achieving a 15% weight reduction (compared to other wraparound and VS models).
- +** **Tri-tab Closure**
Ensures a proper fit.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket
Embroidery	Logo	Sleeve



WRAPAROUND FLEXBACK APRON

MODEL: [AP-WAFB]

Want the ultimate in comfort and protection? The Wraparound Flexback apron combines the best qualities of our flexback apron with the full protection and back support of our wraparound back relief apron. The elastic stretch-back shifts weight from your shoulders to your hips while the tri-tab velcro closure system provides convenient, easy-on easy-off wearability. The Wrap Around Flexback eliminates excess overlapping lead and is up to **15% lighter than traditional wrap styles.**

Core Material	Weight (0.25mm)	Thickness Options
---------------	-----------------	-------------------

Pb	Lead	3400 gsm	0.25, 0.50mm
☝	Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☝	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.