

FULL PROTECTION

## APRON FEATURES

- +** **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- +** **Single Front Panel**  
Eliminates excess overlapping panels achieving a 15% weight reduction (compared to other wraparound and VS models).
- +** **Tri-tab Closure**  
Ensures a proper fit.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Velcro Closures**

### ADDITIONAL CUSTOMIZATION OPTIONS

|                |       |              |
|----------------|-------|--------------|
| Thyroid Collar | Wings | Extra Pocket |
| Embroidery     | Logo  | Sleeve       |

| Core Material | Weight (0.25mm) | Thickness Options |
|---------------|-----------------|-------------------|
|---------------|-----------------|-------------------|

|   |          |                    |
|---|----------|--------------------|
|  Prolite  | 2850 gsm | 0.25, 0.35, 0.50mm |
|  Premier | 2600 gsm | 0.25, 0.35, 0.50mm |
|  Nano    | 2300 gsm | 0.25, 0.35, 0.50mm |



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# WRAPAROUND FLEXBACK MATERNITY APRON

MODEL: [ AP-WAFBMA ]

The Wraparound Flexback Maternity apron combines the best qualities of our flexback apron with the full protection and back support of our wraparound back relief apron. The elastic stretch-back shifts weight from your shoulders to your hips while the tri-tab velcro closure system provides convenient, easy-on easy-off wearability. In addition, the 6" wide frontal velcro allows the wearer to adjust the crossties up or down to account for mom's growing belly. The Wrap Around Flexback eliminates excess overlapping lead and is up to **15% lighter than traditional wrap styles**.

