

FRONT PROTECTION

APRON FEATURES



Backstrap
Velcro backstrap securely fastens shoulder area.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.



Back Relief Belt
Wide, elastic belt provides lumbar support and maximum comfort.

Apron Pocket

Belt Buckle Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Wings
Extra Pocket	Embroidery	Back Mesh
Logo	Sleeve	

Core Material Weight (0.25mm) Thickness Options

Core Material	Weight (0.25mm)	Thickness Options	
Pb	Lead	3400 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm	
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.25, 0.35, 0.50mm	



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

BACK RELIEF APRON

MODEL: [AP-BR]

Save yourself from back pain with Protech's Back Relief apron. This apron features a unique padded belt system constructed with a bonded elasticized panel that will not lose its shape. When engaged, the special back relief belt provides lumbar support while transferring the weight of the apron from your shoulders to your hips for increased comfort and back support throughout the life of the apron. A velcro backstrap is also included to secure the apron around your shoulder area and prevent sliding.

