

FULL PROTECTION



BALANCE™ VEST AND SKIRT APRON

APRON FEATURES

MODEL: [AP-V]_(vest only) [AP-SKIRT]_(skirt only) [AP-VS]_(vest and skirt combo)

- +** Adjustable Buckle Set
Adjustable buckle closures ensure a secure fit.
- +** Non-slip Backing (skirt interior)
Reduces slippage and improves fit
- +** Seam-to-seam Protection
Seam-to-seam, standardized protection in Front (no change in protection at overlap)
- Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket
- Velcro and Buckle Closures
- Hanging Tabs on Skirt

Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Balance configuration **features seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with a premium non-slip fabric in the back.

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar Wings Extra Pocket Suspenders (skirt)
Embroidery Logo Sleeve

Core Material	Weight (0.25mm)	Thickness Options
Lead	3500 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.



Example of a 0.35mm LE Balance Vest Skirt Apron Lead Configuration

