

FULL PROTECTION



BALANCE™ VEST AND SKIRT APRON

APRON FEATURES

MODEL: [AP-V]_(vest only) [AP-SKIRT]_(skirt only) [AP-VS]_(vest and skirt combo)

- ⊕ Adjustable Buckle Set
Adjustable buckle closures ensure a secure fit.
- ⊕ Non-slip Backing (skirt interior)
Reduces slippage and improves fit
- ⊕ Seam-to-seam Protection
Seam-to-seam, standardized protection in Front (no change in protection at overlap)
- Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket
- Velcro and Buckle Closures
- Hanging Tabs on Skirt


Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Balance configuration features **seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with a premium non-slip fabric in the back.

ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Suspenders (skirt)
- Embroidery
- Logo
- Sleeve



Example of a 0.35mm LE Balance Vest Skirt Apron Lead Configuration



Core Material Weight (0.25mm) Thickness Options

Core Material	Weight (0.25mm)	Thickness Options
Pb Lead	3400 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.