

FULL PROTECTION



# BALANCE™ FLEXBACK VEST AND SKIRT APRON

## APRON FEATURES

MODEL: [AP-VFB]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [AP-VSFB]<sub>(vest and skirt combo)</sub>

- + **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- + **Flexback Belt**  
Wraps internally around your back and securely fastens in the front.
- + **Non-slip Backing** (skirt interior)  
Reduces slippage and improves fit.
- + **Seam-to-seam Protection**  
Seam-to-seam, standardized protection in Front (no change in protection at overlap)
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro and Buckle Closures**
- Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar    Wings    Extra Pocket    Suspenders (skirt)  
Embroidery        Logo        Sleeve

The Vest Skirt Flexback shifts 70% of the weight directly to the hips with added support of an adjustable, elastic back relief belt. The Balance configuration **features seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with an anti-slip fabric in the back (interior) to keep your skirt in place.

**Example of a 0.50mm LE Balance Vest Skirt Apron Lead Configuration**

Core Material	Weight (0.25mm)	Thickness Options
Lead	3500 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

