

## FRONT PROTECTION

# APRON FEATURES



### Wide Body Back Panels

The comfortable panels secure your apron and provides lumbar support.

### Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

### Apron Pocket

### Durable Buckle Closure

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Elastic Back
Extra Pocket	Embroidery	Backstrap
Logo	Sleeve	

## Core Material Weight (0.25mm) Thickness Options

Pb	Lead	3500 gsm	0.25, 0.50mm
☐	Prolite	2850 gsm	0.25, 0.35, 0.50mm
IEC	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☐	Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

# BUCKLE APRON

MODEL: [ AP-BA ]

The Buckle Apron comes standard with an easy-release buckle allowing for quick removal of the apron. The belt is adjustable to accommodate different sizes. Criss-cross back straps provide a secure, comfortable fit.

