

FRONT PROTECTION

APRON FEATURES



Wide Body Back Panels
The comfortable panels secure your apron and provides lumbar support.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Durable Buckle Closure

BUCKLE APRON

MODEL: [AP-BA]

The Buckle Apron comes standard with an easy-release buckle allowing for quick removal of the apron. The belt is adjustable to accommodate different sizes. Criss-cross back straps provide a secure, comfortable fit.

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Elastic Back
Extra Pocket	Embroidery	Backstrap
Logo	Sleeve	

Core Material Weight (0.25mm) Thickness Options

Pb	Lead	3400 gsm	0.25, 0.50mm
☑	Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☑	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

