

FULL PROTECTION

## APRON FEATURES

- +** **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- +** **Non-slip Backing (skirt interior)**  
Reduces slippage and improves fit
- +** **Extended Front Panels**  
Front panels are extended to 3/4 overlap coverage where LE protection is doubled
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro and Buckle Closures**
- Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket	Suspenders (skirt)
Embroidery	Logo	Sleeve	

Core Material	Weight (0.25mm)	Thickness Options
---------------	-----------------	-------------------

Core Material	Weight (0.25mm)	Thickness Options
Lead	3400 gsm	0.25, 0.50mm
Prolite	2850 gsm	0.25, 0.35, 0.50mm
Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# ESSENTIAL™ VEST AND SKIRT APRON

MODEL: [AP-VQ]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [AP-VSQ]<sub>(vest and skirt combo)</sub>

Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Essential apron configuration features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection. The skirt comes standard with a premium non-slip fabric in the back. **Essential aprons are only available in standard sizes.**

For a made-to-measure Vest Skirt, see the Balance Vest Skirt Apron.



**Example of a 0.35mm LE Essential Vest Skirt Apron Lead Configuration**

