

FULL PROTECTION



ESSENTIAL™ FLEXBACK VEST AND SKIRT APRON

APRON FEATURES

MODEL: [AP-VFBQ]_(vest only) [AP-SKIRT]_(skirt only) [AP-VSFBQ]_(vest & skirt combo)

- +** **Adjustable Buckle Set**
Adjustable buckle closures ensure a secure fit.
- +** **Flexback Belt**
Wraps internally around your back and securely fastens in the front.
- +** **Non-slip Backing (skirt interior)**
Reduces slippage and improves fit.
- +** **Extended Front Panels**
Front panels are extended to 3/4 overlap coverage where LE protection is doubled
- +** **Shoulder Padding**
Premium foam padding in the shoulders provides extra comfort.
- +** **Apron Pocket**
- +** **Velcro and Buckle Closures**
- +** **Hanging Tabs on Skirt**

The Essential Vest Skirt Flexback shifts 70% of the weight directly to the hips with added support of an adjustable, elastic back relief belt. The Essential apron configuration features **extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection.** The skirt comes standard with an anti-slip fabric in the back (interior) to keep your skirt in place. **Essential aprons are only available in standard sizes.** For a made-to-measure Vest Skirt Flexback, see the Balance Vest Skirt Flexback Apron.

ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Suspenders (skirt)
- Embroidery
- Logo
- Sleeve

Core Material	Weight (0.25mm)	Thickness Options	
Pb	Lead	3400 gsm	0.25, 0.50mm
P	Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
N	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.



Example of a 0.50mm LE Essential Vest Skirt Apron Lead Configuration

