

## FULL PROTECTION

### APRON FEATURES

**⊕ Complete Coverage**  
The front & back panels ensure you are fully protected from beginning to end of your pregnancy.

**⊕ Half Apron Included**  
A Half apron is included for additional protection as you progress through your pregnancy.

**⊕ Tri-tab Closure**  
Ensures a proper fit and easier donning.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**




**Tie & Buckle Closures**

#### ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar   Wings   Extra Pocket  
Embroidery   Logo   Sleeve

The maternity apron features easy ties and a quick release buckle to secure the back and front panels; allowing complete adjustability during pregnancy. The tri-tab velcro closure system also provides convenient, easy-on easy-off wearability. The maternity apron includes an additional 0.50mm equivalent half apron for a total of 1.00mm Pb (when choosing 0.50mm LE) protection over the abdominal area and lower extremities.

#### Core Material   Weight (0.25mm)   Thickness Options

	Prolite	2850 gsm	0.25, 0.35, 0.50mm
	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
	Premier	2600 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

# MATERNITY APRON

MODEL: [ AP-MA ]

