

FULL PROTECTION

## APRON FEATURES



### Complete Coverage

The front & back panels ensure you are fully protected from beginning to end of your pregnancy.



### Half Apron Included

A Half apron is included for additional protection as you progress through your pregnancy.

### Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

### Apron Pocket

### Tie & Buckle Closures

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket
Embroidery	Logo	Sleeve

## Core Material Weight (0.25mm) Thickness Options

Core Material	Weight (0.25mm)	Thickness Options
Prolite	2850 gsm	0.25, 0.35, 0.50mm
Premier	2600 gsm	0.25, 0.35, 0.50mm
Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# MATERNITY APRON

MODEL: [ AP-MA ]

The maternity apron features easy ties and a quick release buckle to secure the back and front panels; allowing complete adjustability during pregnancy. The maternity apron also includes an additional .50mm equivalent half apron for a total of 1.00mm Pb protection over the abdominal area and lower extremities.

