

## FRONT PROTECTION

# APRON FEATURES



### Hanging Straps

Hanging straps on the shoulders allow for easy storage on any Protech apron rack.

### Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.



### Surgical Drop

Quick release velcro closure tabs help maintain sterile field.

### Apron Pocket


### Velcro Closure

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar   Extra Pocket   Embroidery  
Logo

Maintain your sterile field with the quick release Surgical Drop Apron. Easy-release velcro closures allow the surgical drop to fall even when underneath your scrubs. With padded shoulders and a criss-cross back strap design the surgical drop apron provides maximum comfort for all users. Hanging straps are included to allow for easy storage on any Protech apron rack.

Core Material   Weight (0.25mm)   Thickness Options

	Pb	Lead	3500 gsm	0.25, 0.50mm
		Prolite	2850 gsm	0.25, 0.35, 0.50mm
		Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
	PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
		Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

