

## FRONT PROTECTION

# APRON FEATURES

**+** **Hanging Straps**  
Hanging straps on the shoulders allow for easy storage on any Protech apron rack.

**+** **Surgical Drop**  
Quick release velcro closure tabs help maintain sterile field.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Velcro Closure**

MODEL: [ AP-SD ]

Maintain your sterile field with the quick release Surgical Drop Apron. Easy-release velcro closures allow the surgical drop to fall even when underneath your scrubs. With padded shoulders and a criss-cross back strap design the surgical drop apron provides maximum comfort for all users. Hanging straps are included to allow for easy storage on any Protech apron rack.

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar    Extra Pocket    Embroidery  
Logo

Core Material    Weight (0.25mm)    Thickness Options

Pb	Lead	3400 gsm	0.25, 0.50mm
☘	Prolite	2850 gsm	0.25, 0.35, 0.50mm
PR	Premier	2600 gsm	0.25, 0.35, 0.50mm
☘	Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

