

FRONT PROTECTION

APRON FEATURES

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket




Tie Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar Embroidery Wings
Logo Sleeve Backstrap

The Tie Apron is an easy-to-use, affordable option for x-ray protection. Utilizing a tie closure, the Tie Apron's criss crossing straps ensure a secure and comfortable fit. Pediatric sizing available.

Core Material Weight (0.25mm) Thickness Options

 Pb	Lead	3500 gsm	0.25, 0.50mm
	Prolite	2850 gsm	0.25, 0.35, 0.50mm
	Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
	Premier	2600 gsm	0.25, 0.35, 0.50mm
	Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

TIE APRON

MODEL: [AP-TA]

