

FRONT PROTECTION

APRON FEATURES

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Tie Closure

ADDITIONAL CUSTOMIZATION OPTIONS





Thyroid Collar
Logo

Embroidery
Sleeve

Wings
Backstrap

The Tie Apron is an easy-to-use, affordable option for x-ray protection. Utilizing a tie closure, the Tie Apron's criss crossing straps ensure a secure and comfortable fit. Pediatric sizing available.

Core Material Weight (0.25mm) Thickness Options

 Lead	3400 gsm	0.25, 0.50mm
 Prolite	2850 gsm	0.25, 0.35, 0.50mm
 Premier	2600 gsm	0.25, 0.35, 0.50mm
 Nano	2300 gsm	0.25, 0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



**CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.*

TIE APRON

MODEL: [AP-TA]

