

## FRONT PROTECTION

# APRON FEATURES

**⊕ Multi-Panel Design**  
Allows for easy maneuvering and a comfortable fit.

**⊕ Backstrap**  
Velcro backstrap securely fastens shoulder area.






**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Velcro Closure**

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar   Drop Shoulder   Wings  
Extra Pocket   Embroidery   Back Mesh  
Logo   Sleeve

Core Material	Weight (0.25mm)	Thickness Options
 Lead	3500 gsm	0.25, 0.50mm
 Prolite	2850 gsm	0.25, 0.35, 0.50mm
 Prolite IEC	3000 gsm	0.25, 0.35, 0.50mm
 Premier	2600 gsm	0.25, 0.35, 0.50mm
 Nano	2300 gsm	0.35, 0.50mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead and Prolite IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 130kVp diagnostic range.

# UROLOGY APRON

MODEL: [ AP-UA ]

If you sit while performing procedures requiring x-ray, you may be unnecessarily exposing yourself to radiation. The Urology Apron features a three panel system at the bottom of the apron which enables the user to sit comfortably while maintaining maximum coverage from x-ray.

