

FRONT PROTECTION

APRON FEATURES



Multi-Panel Design

Allows for easy maneuvering and a comfortable fit.



Backstrap

Velcro backstrap securely fastens shoulder area.

Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar

Extra Pocket

Logo

Drop Shoulder

Embroidery

Sleeve

Wings

Back Mesh

Core Material

Weight (0.25mm)

Thickness Options

| | | | |
|----|---------|----------|--------------------|
| Pb | Lead | 3400 gsm | 0.25, 0.50mm |
| ☘ | Prolite | 2850 gsm | 0.25, 0.35, 0.50mm |
| PR | Premier | 2600 gsm | 0.25, 0.35, 0.50mm |
| ☘ | Nano | 2300 gsm | 0.25, 0.35, 0.50mm |



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

UROLOGY APRON

MODEL: [AP-UA]

If you sit while performing procedures requiring x-ray, you may be unnecessarily exposing yourself to radiation. The Urology Apron features a three panel system at the bottom of the apron which enables the user to sit comfortably while maintaining maximum coverage from x-ray.

