

RADIATION SAFETY
EYEWEAR
CATALOG

PROTECH
MEDICAL
2024



02

Your eyes are 4 times more sensitive to radiation than other parts of your body. Lead eyewear are intended to be worn by persons who are present in the operating room during medical procedures requiring the use of fluoroscopy. They are intended protect the eyes of the operator from x-ray radiation.





EYEWEAR

QUALITY. SERVICE. INTEGRITY.

Protech Medical is one of the largest US-based leaded eyewear manufacturers in the world. Our licensed opticians handcraft each frame, which are meticulously assembled and inspected with care.

Experience expert craftsmanship with Protech's leaded glasses, goggles, fitovers, faceshields, and prescription eyewear available in a variety of affordable and brand name options.

WHO USES PROTECTIVE GLASSES?

Not too long ago, few doctors used leaded eyewear. With the increase in occurrences of early onset cataracts and other x-ray induced complications, lead glasses are now a staple in fluoroscopic procedures. Protective leaded eyewear is used by roughly **95% of physicians** and staff in the following vocations:



CARDIAC CATH



EP LABS



UROLOGY



INTERVENTIONAL
RADIOLOGY



PAIN
MANAGEMENT



ORTHOPEDIC
SURGERY



PROTECTION SUMMARY

Lead Equivalence - 0.50 or 0.75mm Pb
X-ray Tube Voltage - 50 kV - 150 kV

IEC 61331-1:2014 Regulation (EU) 2016/425
CE 2895

OUR EYEWEAR BRANDS

Ray-Ban

OAKLEY



COSTA

WILEY X

ABACO
POLARIZED

"The typical interventional cardiologist is exposed to 20 and 30 mSv each year. This level of exposure is nearly 10 times greater than the level of exposure received by the entire body."

-- US National Library of Medicine

FRAME SIZING

THE BREAKDOWN

Example frame size: Centrist (55-17-140)*

In the above frame size, the 55 corresponds to "A" (Lens size), the 17 represents "B" (Bridge Size), and 140 is the length of "C" (Temple size). *"B" Size not specified. NOTE: Not all frame sizes are representative of how a frame will fit. It is always best practice to try a frame on when possible. For instance, the Nike Brazen is a "large frame" based on its frame size, yet it fits a majority of men and women's faces.

A

EYE SIZE

Horizontally measured from outermost edges of the lens. Approx 40-62mm.

B

TEMPLE SIZE

Measurement of the distance between two lenses usually measured at the nose bridge.

C

BRIDGE SIZE

Measurement of the temple arms that extend from the front of the frame over your ears. The length of your frame's temples is important to ensure stability and comfort.

D

"B" MEASUREMENT

Measurement of the height of the lenses. *For Bifocal and Progressive lenses, a B measurement of at least 25mm is recommended.



YOUR PERFECT FRAME

WHAT TO LOOK FOR

TEMPLES

The length and curvature of your frame's temples is important to ensure stability and comfort.

TEMPLE GRIPS

Available with several eyewear styles, temple grips can provide comfort and stability for the wearer.

LEAD GLASS

Lead glass lenses provide the main protection against x-ray. The lenses greatly influence weight so it is important to balance maximizing coverage with the need for comfort.

SIDE SHIELDS

Side shields are essential for lateral protection. Some frames come standard with side shields and others are optional

NOSE BRIDGE

A well-designed nose bridge is necessary for a proper fit. Flat nose bridges enable drilling for the addition of loupes.

VENTILATED NOSE PADS

Ventilated nose pads keep the frame comfortably gripped on your face and can reduce fogging.



98 SUPERLITE



| | |
|--------------|-----------------|
| MODEL | RE-98 |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 63g |
| FRAME SIZE | 63-15-120 |
| MATERIAL | Nylon |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Brown |

The 98 Superlite is an ultra light-weight frame designed for comfort. Its high curvature provides enhanced lateral/side protection against x-ray. Rubberized temples provide a secure, comfortable fit.

Features:

- **Rubberized nose pads and temple**
- **Lightweight Nylon Frame**
- **Highly curved lenses for optimal x-ray protection**



99 ULTRALITE



| | |
|--------------|--|
| MODEL | RE-99/99SM/99R |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 67g 59.5g (SM) |
| FRAME SIZE | 61-18-140 56-17-140 (SM) |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 58 (SM) |
| CE | Yes |
| COLORS | Brown (gold), Black, Blue, Silver, Fuschia |

The 99 Ultralites provide a wide area of protection against x-ray. Their high curvature provide adequate lateral protection which eliminates the need for side shields. The 99 comes in 2 sizes and an optional soft rubber nose bridge for additional support and comfort (99R). (Blue only available in regular size; Fuschia available only in small size).

Features:

- **Rubberized temples for added comfort**
- **Optional rubberized nosebridge (99R)**
- **Wide field of view**



9941 ULTRALITE



| | |
|--------------|---|
| MODEL | RE-9941/9941A/994153 |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 64g 60g (994153) |
| FRAME SIZE | 55-22-130 53-22-130 (994153) |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional for regular size) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 58 (994153) |
| CE | Yes |
| COLORS | Silver, Black, Red, Blue, Matte Black |

The 9941 Ultralites come standard with a rubberized nose bridge. They are designed to offer comfort and provide protection for the entire eye area without the need of additional lead side shields. The 9941 is a popular frame and ideal for those who have a high power prescription. Available in 2 sizes and optional adjustable nose bridge for regular sized frame (9941A).

Features:

- **Adjustable nose pads**
- **Rubberized temples for added comfort**
- **Ideal for high prescriptions**



AIRBORNE



| | |
|--------------|-------------------------------------|
| MODEL | RE-AIRB |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 74.7g |
| FRAME SIZE | 57-21-145 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 62 |
| CE | Yes |
| COLORS | Clear/Black |

The Airborne is fully loaded with an Asian fit nose bridge, removable facial cavity, and curved lenses for optimal lateral x-ray protection.

Features:

- **Removable facial cavity foam insert**
- **Rubberized Asian Fit Nose Bridge**
- **Highly curved lenses for optimal x-ray protection**



CENTRIST



| | |
|--------------|--|
| MODEL | RE-CENT |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 88g |
| FRAME SIZE | 55-17-140 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) Extended Side Shields (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Black Crystal |

The Centrist is a rectangular style safety frame made from high quality nylon and has permanent extended side shields for added x-ray protection. The sleek black with clear insert frame has an adjustable nosepiece that provides comfort and a secure fit.

Features:

- **Rubberized & adjustable Asian Fit Nose Bridge**
- **Extended Side Shields**



CIRCUIT



| | |
|--------------|---------------------------------------|
| MODEL | RE-CIRC |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 79g |
| FRAME SIZE | 60-18-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS or SSA (optional) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Silver/Black, Red/Gray, Black/Gray |

The Circuit frame features wrap-around lenses which maximize lateral protection without the use of side shields. Rubberized temples and nose pads provide excellent comfort and grip. The Circuit's modern style combined with its optimized functionality ensure long lasting quality and protection.

Features:

- **Rubberized temples for grip and comfort**
- **Versatile fit for all face types**
- **Wide field of view**



MAKO



| | |
|--------------|--------------------------------|
| MODEL | RE-MAKO |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 69.3g |
| FRAME SIZE | 65-15-130 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Gray/Green |

The Mako frame has highly curved lenses which maximize radiation protection coverage. Integrated side shields equipped with lead strip side shields provide essential lateral x-ray protection.

Features:

- **Rubberized nose pads and temples for grip and comfort**
- **Highly curved frame for optimal x-ray protection**
- **Integrated lateral side shields**



MALIBU



| | |
|--------------|---|
| MODEL | RE-MALIB |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 68.1g |
| FRAME SIZE | 54-20-155 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Matte Blue, Matte Black, Matte Brown, Shiny Black |

The Malibu frame modernizes the classic wayfarer design and is offered in several tasteful color selections.

Features:

- **Affordable and fashion forward**



ONYX



| | |
|--------------|--------------------------|
| MODEL | RE-ONYX |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 71.8g |
| FRAME SIZE | 68-19-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | None |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black, Shiny White |

Deep stylistic curves and wide angle lenses combine style with maximum x-ray protection. The Onyx frame features a ventilated nose bridge for better grip, increased comfort, and reduced fogging.

Features:

- **Highly curved lenses for optimal x-ray protection**
- **Rubberized, ventilated nose bridge**



RAZER



| | |
|--------------|--|
| MODEL | RE-RAZER |
| FRAME TYPE | Wraparound |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 62.5g |
| FRAME SIZE | 58-17-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb SS (attached side shield included) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 59 |
| CE | Yes |
| COLORS | Matte Green, Matte Black, Matte Gray, Matte Blue, Matte Red |

The Razer blends a sleek design into a lightweight and durable frame. Its premium features such as a ventilated nose bridge, rubberized temples, and wide-angle lenses maximize comfort and protection from x-ray. The modern style and use of light-weight materials make the Razer an excellent choice for long-procedures.

Features:

- **Ventilated, rubberized nose bridge**
- **Rubberized Temples for grip and comfort**
- **Removable 0.50mm LE side shield inserts**
- **Lightweight frame**



53 WRAP



| | |
|--------------|---|
| MODEL | RE-53/53A |
| FRAME TYPE | Classic |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 84g |
| FRAME SIZE | 58-16-135 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Glass (included) or 0.75mm Pb Glass (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 62 |
| CE | Yes |
| COLORS | Matte Black, Glossy Black, Tortoise |

The 53 wrap is constructed of durable nylon that exhibits superior flexibility while minimizing potential breakage. Its full saddle bridge provides maximum comfort and evenly distributes weight over the nose. Available with adjustable nose pads (53A).

Features:

- **Optional adjustable nose pads (53A)**
- **Spring Hinge Temples**
- **Built-in, glass side shields**



MODERN METAL

PHASE OUT



| | |
|--------------|-------------------------------------|
| MODEL | RE-002 |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 73g |
| FRAME SIZE | 48-20-135 |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Gold, Light Gunmetal |

The ultra light-weight model 002 features a thin wire frame with optional contoured side shields. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **Optional 0.50mm LE side shields**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



METAL AVIATOR



| | |
|--------------|-------------------------------------|
| MODEL | RE-003 |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 82g |
| FRAME SIZE | 54-17-150 |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 |
| CE | Yes |
| COLORS | Gold |

The Aviator frame is a bold edition to our leaded eyewear lineup. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **Optional 0.50mm LE side shields**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



PETITE METAL



| | |
|--------------|-------------------------------------|
| MODEL | RE-007 |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 63g |
| FRAME SIZE | 52-22-150 |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Gold, Light Gunmetal |

The model 007 is a classic style made of optical steel with bronze electroplate. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **Optional 0.50mm LE side shields**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



CLASSIC METAL



| | |
|--------------|-------------------------------------|
| MODEL | RE-008 |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 64g |
| FRAME SIZE | 54-19-140 |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Gold |

Available in a beautiful light-blue color, this retro style frame is an aesthetically pleasing option for metal wire frames. The classic metal frame come standard with 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **0.50mm LE side shields included**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



553S METALITE



| | |
|--------------|-------------------------------------|
| MODEL | RE-553S |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 62g |
| FRAME SIZE | 50-18-140 |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Silver |

The model 553 is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Its silver color and retro style offer a classic look while non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **0.50mm LE side shields included**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



557S METALITE



| | |
|--------------|-------------------------------------|
| MODEL | RE-557S |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 63g |
| FRAME SIZE | 48-20-135 |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Gold |

The model 557 is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Its gold color and retro style offer a classic look while non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **0.50mm LE side shields included**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



661/662S METALITE



| | |
|--------------|---------------------------------------|
| MODEL | RE-661/662S |
| FRAME TYPE | Metals |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 54g 70g (662S) |
| FRAME SIZE | 52-18-144 52-18-145 (662S) |
| MATERIAL | Optical Steel |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included - 662S) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 62 |
| CE | Yes |
| COLORS | Gold |

The model 662 is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Its gold color and modern shape provide a refreshed look while non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

Features:

- **0.50mm LE side shields included**
- **Adjustable silicone nose pads**
- **Versatile fit for all face types**



90 FITOVER



| | |
|--------------|---|
| MODEL | RE-90 |
| FRAME TYPE | Fitovers |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 112g |
| FRAME SIZE | 66-14-120 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Glass (included) or 0.75mm Pb Glass (optional) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black |

Designed to fit over prescription eyewear, model 90 fitovers are an excellent choice for those requiring protection against x-ray and liquid splashes. The large, clear side shields protect the user from lateral exposure while eliminating the “tunnel” effects of look-alikes.

Features:

- **Wide angle field of view**
- **Large glass side shields help eliminate “tunnel effect”**



FITGUARD ALPHA



| | |
|--------------|--------------------------------|
| MODEL | RE-FGUARDA |
| FRAME TYPE | Fitovers |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 91g |
| FRAME SIZE | 69-17-145 |
| MATERIAL | Rubberized Matte Black Finish |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black |

The Fitguard Alpha fitover is a versatile, high quality nylon frame with blacked-out side shields and a streamlined look. Featuring Flexfit technology, the fitover's temples can be easily adjusted for a secure, comfortable fit. The Alpha (XL) model is ideal for covering wide, rectangular shaped prescription eyewear. With balanced weight distribution and the most coverage, the Alpha is our most popular fitover model.

Features:

- **Adjustable temples**
- **Highly durable frame**
- **Ideal for wide prescription eyewear**



FITGUARD BETA



| | |
|--------------|--------------------------------|
| MODEL | RE-FGUARDB |
| FRAME TYPE | Fitovers |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 108g |
| FRAME SIZE | 64-18-150 |
| MATERIAL | Rubberized Matte Black Finish |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black |

The Fitguard Beta fitover is a versatile, high quality nylon frame with blacked-out side shields and a streamlined look. Featuring Flexfit technology, the fitover's temples can be benched (adjusted) easily for a secure, comfortable fit. The Beta (LG) model is ideal for covering taller, more square shaped prescription eyewear.

Features:

- **Adjustable temples**
- **Highly durable frame**
- **Ideal for fitting over eyewear with square or rounded lenses**



FULL FACE SHIELD



| | |
|--------------|------------------------------------|
| MODEL | 450FPM |
| FRAME TYPE | Face Shields |
| FRONT LENS | 0.10mm Pb Acrylic |
| WEIGHT | 350g |
| FRAME SIZE | Adjustable 20" - 25" circumference |
| MATERIAL | Acrylic, Plastic |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Black/Clear |

The Full Face Shield uses a lead impregnated protective acrylic shield suspended from a dual adjustable headpiece. The Full face shield provides maximum x-ray protection due to its larger lead impregnated shield. Because of the multi-adjustable headpiece, you can angle the mask to protect from radiation in many different positions as well as adjust the size to comfortably fit your head. *Mask should be worn with leaded eyewear protection.



Made in the USA
from imported parts

PANORAMIC FACE SHIELD



| | |
|--------------|------------------------------------|
| MODEL | 400PSM |
| FRAME TYPE | Face Shields |
| FRONT LENS | 0.10mm Pb Acrylic |
| WEIGHT | 425g |
| FRAME SIZE | Adjustable 20" - 25" circumference |
| MATERIAL | Acrylic, Plastic |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Black/Clear |

The Panoramic Face Shield utilizes a lead impregnated acrylic shield suspended from a dual adjustable headpiece. This lead face shield provides additional x-ray protection around the sides of the head while sacrificing some front protection. Because of the multi-adjustable headpiece, you can angle the mask to protect from radiation in many different positions as well as adjust the size to comfortably fit your head. *Mask should be worn with leaded eyewear protection.



Made in the USA
from imported parts

70 ASTROFLEX



| | |
|--------------|-------------------------------------|
| MODEL | RE-70S |
| FRAME TYPE | Goggles |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 88g |
| FRAME SIZE | 54-19-Adjustable |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 62 |
| CE | Yes |
| COLORS | Clear/Blue |

Impact resistance and flexibility come together in this high quality goggle. The frame features molded side shielding and is fitted to the brow for increased splash protection. The temple length and frame angle are adjustable for a uniform custom fit. Model 70S comes with leaded side shields for necessary lateral x-ray coverage.

Features:

- Adjustable temples
- Wide field of view



COMET ECONOMY



| | |
|--------------|--|
| MODEL | RE-COM |
| FRAME TYPE | Goggles |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 65g |
| FRAME SIZE | 50-19-117 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) Extended Side shield Seal (optional) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Orange/Black, Yellow/Black, Black |

The Comet Economy is an affordable option providing excellent coverage against x-ray. The side shield configuration allows for lead shielding to be placed on the frame exterior and wrap around the entire lens area, providing a complete seal from scatter radiation.

Features:

- **Rubberized nose-fit bridge and temples for grip and comfort**
- **Foam top bar for added comfort**
- **Optional Extended side shields**
- **Wide field of view**



COMET



| | |
|--------------|--|
| MODEL | RE-COMET |
| FRAME TYPE | Goggles |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 78.3g |
| FRAME SIZE | 59-19-114 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) Extended Side shield Seal (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 |
| CE | Yes |
| COLORS | Blue, Gray/Orange, Pink |

A sleek and versatile relative to the Proguard LeadR's frame, the Comet provides excellent coverage against x-ray. Its unique side shield configuration allows for lead shielding to be placed on the frame exterior and wrap around the entire lens area, providing a complete seal from scatter radiation.

Features:

- **Rubberized nose-fit bridge and temples for grip and comfort**
- **Foam top bar for added comfort**
- **Optional Extended side shields**
- **Wide field of view**



LEAD'R



| | |
|--------------|--|
| MODEL | RE-OG220 |
| FRAME TYPE | Goggles |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 76.7g |
| FRAME SIZE | 58-15-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (included) |
| | Extended Side Shields (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 66 |
| CE | Yes |
| COLORS | Black, Gray/Pink, Green/Yellow, White/Purple |

The LeadR's feature a soft, ventilated nose bridge with rubberized temples for premium comfort. Each LeadR's frame comes with the added protection of 0.50mm leaded side shields and the ability to extend that side protection into an additional panel.

Features:

- **Wide field of view**
- **Rubberized, ventilated nose bridge**
- **Rubberized, curved temples**
- **Foam top bar for added comfort**
- **Optional Extended side shields**



X-GUARD



| | |
|--------------|-------------------------------------|
| MODEL | RE-XGUARD |
| FRAME TYPE | Goggles |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 81.8g |
| FRAME SIZE | 58-18-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips (included) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 62 |
| CE | Yes |
| COLORS | Black/Woodgrain, Matte Black |

Tired of the basic goggle selections for radiation protection? So are we! The X Guard goggle features a premium, clear insert, lead strip side shields for lateral coverage (shown in black model below), and tasteful color options.

Features:

- **Wide field of view**
- **Removable side shield insert**



ABACO DOCKSIDE



| | |
|--------------|-------------------------------------|
| MODEL | AB-DOCK |
| FRAME TYPE | Brand - Abaco |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 72g |
| FRAME SIZE | 58-19-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 60 |
| CE | Yes |
| COLORS | Black, Crystal Clear |

Dockside glasses have a bold square style in an impact-resistant polycarbonate frame featuring premium metal hinges. These fashion forward frames are both durable and lightweight.

Features:

- **Lightweight frame**
- **Versatile fit for all face types (unisex)**



ABACO JESSE



| | |
|--------------|-------------------------------------|
| MODEL | AB-JES |
| FRAME TYPE | Brand - Abaco |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 68g |
| FRAME SIZE | 56-18-142 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 60 |
| CE | Yes |
| COLORS | Matte Black |

These unisex, universally flattering eyewear are made of soft, lightweight polycarbonate. They feature characteristic rivets at the temples and extremely flexible spring hinges for a snug, comfortable fit.

Features:

- **Lightweight frame**
- **Versatile fit for all face types (unisex)**



ABACO KAI



| | |
|--------------|-------------------------------------|
| MODEL | AB-KAI |
| FRAME TYPE | Brand - Abaco |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 65g |
| FRAME SIZE | 56-18-135 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 60 |
| CE | Yes |
| COLORS | Crystal Clear, Matte Black |

Made of shatter-resistant, lightweight polycarbonate, the frames feature premium metal hinges and are molded around the ears to provide a comfortable fit. These unisex sunglasses complement a variety of facial shapes, and they fit virtually all faces from small to large.

Features:

- **Lightweight frame**
- **Versatile fit for all face types (unisex)**



ABACO WAIKIKI



| | |
|--------------|-------------------------------------|
| MODEL | AB-WAI |
| FRAME TYPE | Brand - Abaco |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 68g |
| FRAME SIZE | 54-21-140 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 60 |
| CE | Yes |
| COLORS | Tortoise, Gloss Black |

This revived version of an iconic classic features a lightweight polycarbonate frame to provide durability and a comfortable fit. These sunglasses feature characteristic rivet details on the front, premium metal hinges, and they are molded around the ears for added comfort. These eyewear look good on virtually everyone.

Features:

- **Lightweight frame**
- **Versatile fit for all face types (unisex)**



COSTA ALETA



| | |
|--------------|-------------------------------------|
| MODEL | CO-ALE |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 54-19-136 |
| MATERIAL | Bioplastics |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Wahoo, Taupe Crystal |

Offering performance in a stylish yet functional design, Aleta uses lightweight Bio-Resin and double-injected Hydrolite™ for a better grip. With micro shields and hooding, excess light stays out while fog stays down with vented nose pads.

Features:

- **Rubberized nose pads**
- **Lightweight and durable plant-based frame (bioplastics)**



COSTA ARANSAS



| | |
|--------------|-------------------------------------|
| MODEL | CO-ARA |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 58-17-145 |
| MATERIAL | Acetate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Ocean Tortoise |

The Aransas features an acetate frame with a classic square lens shape and crafted design elements. This versatile frame fits most face shapes and offers medium coverage.

Features:

- **High-fit Bridge**
- **Universal Fit (for most face shapes)**



COSTA BLACKFIN



| | |
|--------------|-------------------------------------|
| MODEL | CO-BFIN |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 62-17-120 |
| MATERIAL | Bioplastics |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Matte Black |

Flexible and boasting a Hydrolite™ co-injected lining, the Costa Blackfin lead glasses are made for performance and are a durable option for any doctor.

Features:

- **High-fit Bridge**
- **Rubberized Nosepads**
- **Lightweight and durable plant-based frame (bioplastics)**



COSTA BIMINI ROAD 330



| | |
|--------------|-------------------------------------|
| MODEL | CO-BRD330 |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 55-16-145 |
| MATERIAL | Metal |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Black, Pacific Blue, Gunmetal |

The Bimini Road features a smaller lens designed to fit those with a narrow head. Classic style designed in titanium and Monel®, these durable, lightweight frames are corrosion-resistant and a reliable choice for anyone.

Features:

- **Adjustable nose pads**
- **Narrow fit for smaller faces**
- **Temporary Removable Side Shields (optional)**
- **Lightweight, durable, corrosion resistant frame**



COSTA BIMINI ROAD 320



| | |
|--------------|-------------------------------------|
| MODEL | CO-BRD320 |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 53-17-145 |
| MATERIAL | Metal |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Black, Pacific Blue, Gunmetal |

The Bimini Road features a smaller lens designed to fit those with a narrow head. Classic style designed in titanium and Monel®, these durable, lightweight frames are corrosion-resistant and a reliable choice for anyone.

Features:

- **Adjustable nose pads**
- **Narrow fit for smaller faces**
- **Temporary Removable Side Shields (optional)**
- **Lightweight, durable, corrosion resistant frame**



COSTA FERNANDINA



| | |
|--------------|-------------------------------------|
| MODEL | CO-FER |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 57-15-132 |
| MATERIAL | Metal, Acetate |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Brushed Gold |

The Fernandina is a small metal frame with rounded lenses and acetate temples for added color. This unisex frame for smaller faces has premium 5-barrel hinges and adjustable nosepads for added comfort.

Features:

- **High-fit Bridge**
- **Narrow fit for smaller faces**
- **Adjustable Nosepads**
- **Temporary Removable Side Shields (optional)**





COSTA OCEAN RIDGE 400

| | |
|--------------|-------------------------------------|
| MODEL | CO-OCR4 |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 57-16-140 |
| MATERIAL | Bioplastics |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Black, Pacific Blue, Fog Gray |

The Ocean Ridge features a lightweight frame made from Bio-Resin™ nylon. The material is strong, allowing the frames to hold their shape in heat and cold weather. This versatile frame fits most face shapes.

Features:

- **Universal Fit (for most face shapes)**
- **Rubberized Nosepads**
- **Temporary Removable Side Shields (optional)**
- **Lightweight and durable plant-based frame (bioplastics)**



COSTA SEADRIFT

| | |
|--------------|---------------------|
| MODEL | CO-SEA |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 60-15-123 |
| MATERIAL | Bioplastics |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | No |
| COLORS | Shiny Tortoise Fade |

This sleek, Bio-Resin frame sports Costa's Hydrolite™ rubber in the temple arms and nose pads. Providing optimal coverage and wrap, the Seadrift is there when you need to perform.

Features:

- **High-fit Bridge**
- **Rubberized Nosepads and temples**
- **Lightweight and durable plant-based frame (bioplastics)**



COSTA SULLIVAN



| | |
|--------------|--|
| MODEL | CO-SUL |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 53-19-140 |
| MATERIAL | Acetate |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Tortoise, Black Kelp, Matte Gray Crystal |

Sullivan was designed with a slow-down, kick-back kind of style along with its water-inspired acetate designs. This frame fits all face types and offers medium coverage.

Features:

- **High bridge fit**
- **Temporary Removable Side Shields (optional)**
- **Versatile fit for all face types (unisex)**



COSTA TUNA ALLEY

| | |
|--------------|-----------------|
| MODEL | CO-TALY |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 61.5-11-125 |
| MATERIAL | Nylon |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | No |
| COLORS | Blackout |

With integral hinges, a patented vent system, and durable nylon construction, these men's performance glasses won't be missing much. The Tuna Alley is a large frame offering maximum coverage and wrap.

Features:

- **Rubberized nose pads**
- **High bridge and wide fit**
- **Durable Nylon Frame**
- **Vent System**
- **Maximum coverage**



COSTA TYBEE



| | |
|--------------|--|
| MODEL | CO-TYB |
| FRAME TYPE | Brand - Costa |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 55-20-140 |
| MATERIAL | Acetate |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | - |
| CE | No |
| COLORS | Tortoise, Matte Black, Ocean Currents, Shiny Light Gray Crystal |

Boasting ocean-inspired colors and a casual style, the Tybee features a standard fit with medium coverage and wrap.

Features:

- **High-fit Bridge**
- **Temporary Removable Side Shields (optional)**



NIKE 7118



| | |
|--------------|---------------------------------------|
| MODEL | N-7118 |
| FRAME TYPE | Brand - Nike |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 71.5g |
| FRAME SIZE | 55-14-140 |
| MATERIAL | Bioplastics |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 55 |
| CE | No |
| COLORS | Matte Midnight Navy, Matte Black/Blue |

The Nike 7118 lead glasses are a retro-chic, gender-neutral frame made of lightweight, bio-plastic. Featuring rubberized nose pads and temple tips, the 7118 protective eyewear are a comfortable "all-day-wear" option.

Features:

- **Rubberized nose pads and temple tips**
- **Temporary Removable Side Shields (optional)**
- **Lightweight and durable plant-based frame (bioplastics)**



NIKE BRAZEN



| | |
|--------------|-------------------------------------|
| MODEL | N-BR |
| FRAME TYPE | Brand - Nike |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 75g |
| FRAME SIZE | 60-16-140 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 66 (Min. Single Vision PD: 64) |
| CE | Yes |
| COLORS | Matte Black with Blue, Shiny Black |

The Nike Brazen features an advanced ventilated nose bridge for grip, comfort, and reduced fogging. Its lightweight nylon design combines streamlined style with enhanced comfort. The premium Nike Brazen frame offers versatility with its superior coverage and protection and fits a wide range of people.

Features:

- **Rubberized, ventilated nose bridge**
- **Lightweight frame**
- **Superior coverage for optimal x-ray protection**



NIKE BRAZEN FURY



| | |
|--------------|---|
| MODEL | N-BRF |
| FRAME TYPE | Brand - Nike |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | - |
| FRAME SIZE | 60-17-130 |
| MATERIAL | Bioplastics |
| SIDE SHIELDS | 0.50mm Pb Strips SS or SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 |
| CE | No |
| COLORS | Matte Black with Green, Matte Gray/Blue |

The Nike Brazen Fury features rubberized nose pads and temple tips for added grip and comfort. With super-comfortable features in a lightweight frame, they're an ideal pair for active and everyday wear.

Features:

- **Rubberized nose pads and temple tips**
- **Lightweight and durable plant-based frame (bioplastics)**



NIKE MAVERICK



| | |
|--------------|-------------------------------------|
| MODEL | N-MAV |
| FRAME TYPE | Brand - Nike |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 70.3g |
| FRAME SIZE | 59-15-145 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 |
| CE | No |
| COLORS | Matte Black |

Utilizing Nike's unique nose bridge design, the Maverick frame reduces fogging and maximizes comfort during long fluoroscopic procedures.

Features:

- **Unique nose bridge design for fog reduction**
- **Lightweight frame**



NIKE PREMIER

| | |
|--------------|---|
| MODEL | N-PRMR |
| FRAME TYPE | Brand - Nike |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 67.8g |
| FRAME SIZE | 60-13-135 |
| MATERIAL | Nylon |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 64 |
| CE | No |
| COLORS | Matte Tortoise, Matte Anthracite, Matte Black |

The Nike Premier is reborn with a diverse selection of colors and tasteful design elements. The Premier's large rectangular lenses provide adequate x-ray protection and coverage.

Features:

- **Large, rectangular lenses**
- **Lightweight frame**



NIKE RABID



| | |
|--------------|---------------------------------|
| MODEL | N-RABID |
| FRAME TYPE | Brand - Nike |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 63g |
| FRAME SIZE | 63-14-145 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | None |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black/White |

The Rabid is a high-curvature sports frame providing superior wraparound coverage and protection. Its lightweight nylon design and wide lenses provides clarity at all angles during any procedure. Secure Wrap Technology features grip zones throughout the frame for a better overall fit, comfort, and stability.

Features:

- **High curvature frame for wraparound coverage and optimal x-ray protection**
- **Lightweight frame**
- **Superior grip for stability and comfort**



LIBERTY MAXX 30



| | |
|--------------|--|
| MODEL | LBT-MAX30S (small) LBT-MAX30 (reg.) |
| FRAME TYPE | Brand - Liberty |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 60g (small) 69g (reg.) |
| FRAME SIZE | 53-17-130 (small) 55-20-130 (reg.) |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 56 |
| CE | Yes |
| COLORS | Shiny Black, Shiny Blue, Shiny Red, Shiny Silver |

The Liberty Maxx's hingeless temples are ultra-sturdy while the head strap keeps these glasses fitting securely during long procedures. Soft padding inside the temples and bridge ensures maximum comfort and the best possible fit.

Features:

- **Available in 2 sizes**
- **Hingeless temples for a secure fit**
- **Can accommodate wide range of prescriptions**
- **Versatile fit for all face types (unisex)**
- **Rubberized nose bridge and temples for added comfort and fit**



LIBERTY MAXX 31



| | |
|--------------|---------------------------------------|
| MODEL | LBT-MAX31S (small) LBT-MAX31 (reg.) |
| FRAME TYPE | Brand - Liberty |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 60g (small) 64g (reg.) |
| FRAME SIZE | 53-17-130 (small) 55-20-130 (reg.) |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 56 |
| CE | Yes |
| COLORS | Shiny Blue, Shiny Black, Shiny Silver |

The Liberty Maxx's head strap keeps these glasses fitting securely during long procedures. Soft padding inside the temples and bridge ensures maximum comfort and the best possible fit.

Features:

- Available in 2 sizes
- Shortened temples with a fitted head strap
- Can accommodate wide range of prescriptions
- Versatile fit for all face types (unisex)
- Rubberized nose bridge for comfort



OAKLEY CROSSLINK ZERO



| | |
|--------------|-------------------------------------|
| MODEL | OA-CROSS |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 61.6g |
| FRAME SIZE | 58-16-138 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 56 |
| CE | Yes |
| COLORS | Matte Black, Polished Gray Shadow |

The Crosslink Zero has fixed temples with a three-Point Fit to ensure lenses are aligned for visual correction. Wire cores hidden inside the temples of this lightweight frame offer adjustability for a perfect fit.

Features:

- Rubberized nose pads and ear socks for added comfort
- Lightweight frame



OAKLEY FIVE SQUARED



| | |
|--------------|----------------------------|
| MODEL | OA-FIVESQ |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 71g |
| FRAME SIZE | 54-20-135 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Polished Black, Gray Smoke |

Oakley Five Squared leaded glasses uses Condensed Cranial Geometry for a tapered architecture that feels as good as it looks and it's a look that soaks sophistication in adrenalin. Dimensional reliefs, metal icons and sculpturally integrated hinges take the whole idea of streetsmooth style and amp it with attitude.

Features:

- Lightweight frame



OAKLEY HOLBROOK



| | |
|--------------|--|
| MODEL | OA-HOL |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 66.5g |
| FRAME SIZE | 56-18-135 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 60 |
| CE | No |
| COLORS | Polished Clear, Satin Black, Matte Brown Tortoise |

The Oakley Holbrook frame combines traditional design with modern style and flare.

Features:

- Lightweight frame



OAKLEY JUPITER SQUARED



| | |
|--------------|-------------------------------------|
| MODEL | OA-JUP |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 74.1g |
| FRAME SIZE | 56-18-135 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 66 |
| CE | Yes |
| COLORS | Matte Black - USA Flag |

A customer favorite, the sporty Jupiter Squared is attitude with refinement. The handles of switchblade knives inspired the earstems and their metal bands and rivets. The Matte Black frames of the Flag Collection are highlighted with American flag emblems.

Features:

- **Rubberized ear socks for added grip and comfort**
- **American flag on Oakley "O" Logo**



OAKLEY LITEBEAM



| | |
|--------------|-------------------------------------|
| MODEL | OA-LITE |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 59.6g |
| FRAME SIZE | 55-18-140 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 60 |
| CE | No |
| COLORS | Satin Black |

Smooth details and streamlined design provide a timeless look with the Oakley Litebeam. Composed of lightweight materials, the Litebeam features adjustable temples and comes with several nosebridge sizes to ensure maximum comfort.

Features:

- **Lightweight frame**
- **Adjustable temples for a perfect fit**



OAKLEY STRAIGHTLINK

PHASE OUT



| | |
|--------------|-------------------|
| MODEL | OA-STRGHT |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 77.2g |
| FRAME SIZE | 61-19-130 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black |

Unobtainium® earsocks and nose pads increase grip with perspiration for a no slip grip. Three-Point Fit holds lenses in precise optical alignment while eliminating pressure points. O Matter™ stress-resistant frame material is both lightweight and durable for all-day comfort and protection.

Features:

- **Lightweight frame**
- **Optimised weight distribution and alignment**
- **Rubberized ear socks and nose pads for comfort and grip**



OAKLEY TRILLBE X



| | |
|--------------|-------------------------------------|
| MODEL | OA-TRILL |
| FRAME TYPE | Brand - Oakley |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 64g |
| FRAME SIZE | 52-18-140 |
| MATERIAL | Durable O Matter™ |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 58 |
| CE | Yes |
| COLORS | Polished Clear |

Designed for small to medium faces, this comfortably lightweight addition to the Trillbe™ Collection has a softly rounded front with a sculptured keyhole bridge, and its stem architecture takes vintage style into modern territory.

Features:

- **Lightweight frame**
- **Designed for small to medium faces**



RAY BAN 2016 DADDY-O



| | |
|--------------|-------------------------------------|
| MODEL | RB-2016 |
| FRAME TYPE | Brand - Ray Ban |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 61g |
| FRAME SIZE | 59-17-125 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Strips SS (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 65 |
| CE | Yes |
| COLORS | Matte Black |

Go 100% classic Ray-Ban with the laid back, classic style of Daddy-O.

Features:

- **Lightweight frame**
- **Excellent coverage for maximum x-ray protection**



RAY BAN 4105



| | |
|--------------|-------------------------------------|
| MODEL | RB-4105 |
| FRAME TYPE | Brand - Ray Ban |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 62g |
| FRAME SIZE | 50-23-140 |
| MATERIAL | Nylon |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 65 |
| CE | Yes |
| COLORS | Matte Black |

The Ray-Ban Classic Wayfarer look with a twist: the glasses are foldable. Modeled after Ray-Ban's iconic design, the 4105 was reconstructed into a compact, folding, light weight frame.

Features:

- **Impact resistant, durable nylon**
- **Compact folding style for ease of portability**



RAY BAN 4165F JUSTIN



| | |
|--------------|-------------------------------------|
| MODEL | RB-4165 |
| FRAME TYPE | Brand - Ray Ban |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 73g |
| FRAME SIZE | 55-18-145 |
| MATERIAL | Nylon |
| SIDE SHIELDS | N/A |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 65 |
| CE | Yes |
| COLORS | Matte Black |

Bold, Fun and Unforgettable, the 4165F glasses may just be one of the coolest looks in the Ray-Ban collection. Inspired by the iconic Wayfarer sunglasses, this style has been redesigned with slightly larger rectangle lenses, fresh new frame colors, and a rubberized frame finish to really give it that extra edge.

Features:

- Rubberized frame finish for no-slip grip and comfort
- Durable nylon



RAY BAN 4181



| | |
|--------------|-------------------------------------|
| MODEL | RB-4181 |
| FRAME TYPE | Brand - Ray Ban |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 71.4g |
| FRAME SIZE | 57-16-145 |
| MATERIAL | Injected Plastic |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | Single Vision, Bifocal, Progressive |
| MIN PRG PD | 65 |
| CE | No |
| COLORS | Tortoise |

The stylistically distinct and unassuming Ray Ban model 4181 serves as a timeless option for radiation protection with class.

Features:

- Versatile fit for all face types (unisex)



RAY BAN 7047

| | |
|--------------|---------------------------------|
| MODEL | RB-7047 |
| FRAME TYPE | Brand - Ray Ban |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 63g |
| FRAME SIZE | 54-17-140 |
| MATERIAL | Injected Plastic |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | No |
| COLORS | Black, Blue, Brown |

A sharper more self assured Ray Ban style, the model 7047 is full of confidence and style. Composed of lightweight materials, the 7047 provides comfortable, long lasting protection against x-ray.

Features:

- **Wide fit temples for a comfortable fit**



WILEY-X AIRRAGE

| | |
|--------------|-----------------|
| MODEL | WX-AIRR |
| FRAME TYPE | Brand - Wiley-X |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 67.5g |
| FRAME SIZE | 61-18-124 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | No |
| COLORS | Matte Black |

The Wiley X Airrage frame features a facial foam cavity seal and high curved lenses for maximum comfort and protection against x-ray.

Features:

- **Facial foam cavity for added comfort**
- **Highly curved lenses for optimal x-ray protection**
- **Durable frame**



WILEY-X BOSS



| | |
|--------------|-----------------|
| MODEL | WX-BOSS |
| FRAME TYPE | Brand - Wiley-X |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 84.5g |
| FRAME SIZE | 68-18-125 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black |

Featuring highly curved lenses, the Wiley X Boss provides excellent x-ray protection.

Features:

- **Removable facial foam cavity seal for added comfort**
- **Highly curved lenses for optimal x-ray protection**
- **Durable frame**
- **Maximum splash protection**



Facial Cavity Seal



WILEY-X BRICK



| | |
|--------------|-----------------|
| MODEL | WX-BRICK |
| FRAME TYPE | Brand - Wiley-X |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 74g |
| FRAME SIZE | 63-18-120 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black |

The Wiley X Brick is part of a long line of military inspired frames with sleek, modern designs. The Wiley X Brick's highly curved 8 base lenses provide excellent coverage against x-ray and feature a removable facial cavity seal to protect against debris and splash exposure.

Features:

- **Removable facial foam cavity seal for added comfort**
- **Highly curved lenses for optimal x-ray protection**
- **Durable frame**
- **Rubberized, flexible temple arms**



Facial Cavity Seal



WILEY-X CONTOUR



| | |
|--------------|---------------------------------------|
| MODEL | WX-CON |
| FRAME TYPE | Brand - Wiley-X |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 62.3g |
| FRAME SIZE | 54-17-140 |
| MATERIAL | Nylon |
| SIDE SHIELDS | 0.50mm Pb Inserts SS (included) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Matte Black, Gloss Black/Brown Stripe |

With the move to branded lead glasses, the Wiley X Contour is the whole package. Its retro style and modern color palette deliver on appeal, while optional premium side shield inserts are also available for maximum coverage against x-ray.

Features:

- **Removable 0.50mm side shield inserts included**
- **Durable, lightweight frame**



WILEY-X P-17



| | |
|--------------|---------------------------------|
| MODEL | WX-P17 |
| FRAME TYPE | Brand - Wiley-X |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 68g |
| FRAME SIZE | 60-18-125 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | 0.50mm Pb Strips SSA (optional) |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Gloss Black, Matte Black, Brown |

The Wiley X P-17 is a frame inspired by the military. Its sleek design and lightweight comfort provide unparalleled performance. The Wiley X P-17 is ideal for medium to large facial structures and comes standard with high curvature 0.75mm LE lenses.

Features:

- **Ideal for medium to large face types**
- **Highly curved lenses for optimal x-ray protection**
- **Durable frame**



TUGGA BABY



| | |
|--------------|-------------------------------------|
| MODEL | RE-BABY |
| FRAME TYPE | Goggles |
| FRONT LENS | 0.75mm Pb Glass |
| WEIGHT | 42g |
| FRAME SIZE | 46-18 |
| MATERIAL | Polycarbonate |
| SIDE SHIELDS | N/A |
| RX OPTIONS | N/A |
| MIN PRG PD | N/A |
| CE | Yes |
| COLORS | Shiny Pink, Shiny Blue, Shiny Green |

Protech's Tugga Baby pediatric glasses provide essential protection during medical fluoroscopy and feature a universal fit with its adjustable headband.

Features:

- **Adjustable headband for universal fit**



EYEWEAR ADDONS

LATERAL SIDE SHIELDS

Many of our frames come with side shields for maximum protection. There are three types of side shields - glass, side shield attachments (SSA) and vinyl strips to existing side shields (SS). **For all side shields types, the standard protection is 0.50mm Pb.** Some frames are able to accommodate a 0.75mm Pb glass, such as the 53 Wrap shown here (Fig A.).



Fig A - Side shield example - 0.75mm Pb Glass



Side shield Attachment (SSA) example



Side shield (SS) example - 0.50mm Pb strip

AR COATING

An anti-reflective coating eliminates the reflection produced by light on eyewear lenses. **Once an AR coat has been applied, 99.5% of available light will be able to pass through your lenses (compared to 85%) which enhances visual acuity and reduces eye strain.** During procedures, an anti-reflective coating is particularly important to reduce the reflections and light disturbances caused by operational and lighting equipment.

ACCESSORIES



All eyewear purchases come with a hard case, pouch and cleaning cloth.



RE-CLEAN

RE-CLOTH

RE-CASE



RE-STRAP



RE-STRAP2

| ITEM CODE | DESCRIPTION |
|-----------|--|
| RE-CASE | Honeycomb Hard Case |
| RE-STRAP | Simple Square Slider Eyewear Strap |
| RE-STRAP2 | Lobster Claw Eyewear Strap |
| RE-CLEAN | Anti-fog Eyewear Lens Cleaner 1 fl. oz |
| RE-CLOTH | Soft Eyewear Cloth Wipe |

| | | | | | | | | Prescription (Rx) Options | | | | |
|----------------------|-----------|------------|------------|--------------|----|---------------|------------|---------------------------|----|-------------|----------------|----------------------|
| Model | Size | Base Curve | Weight (g) | Side Shields | CE | Material | Frame Type | Plano | Rx | Min. PRO PD | Powers | Included Accessories |
| Protech Frames | | | | | | | | | | | | |
| Fitovers | | | | | | | | | | | | |
| 90 Fitover | 66-14-120 | 4BC 80 | 112 | Incl | Y | Nylon | Fitover | ✓ | - | - | - | OEM: C, ST, CL, P |
| FitGuard Alpha | 69-17-145 | 4BC 80 | 91 | Incl | Y | Nylon | Fitover | ✓ | - | - | - | OEM: C, ST, CL, P |
| FitGuard Beta | 64-18-150 | 4BC 80 | 108 | Incl | Y | Nylon | Fitover | ✓ | - | - | - | OEM: C, ST, CL, P |
| Goggles | | | | | | | | | | | | |
| 70s Astroflex | 54-19-ADJ | 4BC 65 | 88 | Incl | Y | Polycarbonate | Goggle | ✓ | ✓ | 62 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Comet | 59-19-114 | 6BC 65 | 78.3 | Incl | Y | Polycarbonate | Goggle | ✓ | ✓ | 64 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Comet Economy | 50-19-117 | 6BC 65 | 65 | Incl | Y | Polycarbonate | Goggle | ✓ | - | - | - | OEM: C, ST, CL, P |
| LeadR | 58-15-135 | 6BC 65 | 76.7 | Incl | Y | Polycarbonate | Goggle | ✓ | ✓ | 66 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| X Guard | 58-18-135 | 6BC 65 | 81.8 | Incl | Y | Polycarbonate | Goggle | ✓ | ✓ | 62 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Classic | | | | | | | | | | | | |
| 53 Wrap | 58-16-135 | 4BC 65 | 84 | Incl | Y | Nylon | Glasses | ✓ | ✓ | 62 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Malibu Wayfarer | 54-20-155 | 4BC 65 | 68.1 | N/A | Y | Polycarbonate | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Classic Metals | | | | | | | | | | | | |
| 553s Metalite | 50-18-140 | 4BC 65 | 62 | Incl | Y | Optical Steel | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 557s Metalite | 48-20-135 | 4BC 65 | 63 | Incl | Y | Optical Steel | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 661 Metalite | 52-18-144 | 4BC 65 | 54 | N/A | Y | Optical Steel | Glasses | ✓ | ✓ | 62 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 662s Metalite | 52-18-145 | 4BC 65 | 70 | Incl | Y | Optical Steel | Glasses | ✓ | ✓ | 62 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Classic Metal | 54-19-140 | 4 - 6 | 64 | Incl | Y | Optical Steel | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Metal Aviators | 54-17-150 | 4 - 6 | 82 | SS | Y | Optical Steel | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Modern Metal | 48-20-135 | 4 - 6 | 73 | SS | Y | Optical Steel | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Petite Metal | 52-22-150 | 4 - 6 | 63 | SS | Y | Optical Steel | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Wrap-Around | | | | | | | | | | | | |
| 98 | 63-15-120 | 8BC 73 | 63 | N/A | Y | Nylon | Glasses | - | - | - | - | OEM: C, ST, CL, P |
| 99 Ultralite | 61-18-140 | 6BC 65 | 67 | SS | Y | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 99 Ultralite Small | 56-17-140 | 6BC 65 | 59.5 | N/A | Y | Nylon | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 99 Ultralite Bridge | 61-18-140 | 6BC 65 | 67 | SS | Y | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 9941 Ultralite | 55-22-130 | 6BC 65 | 64 | SS | Y | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| 9941-53 Ultralite SM | 53-22-130 | 6BC 65 | 58 | N/A | Y | Nylon | Glasses | ✓ | ✓ | 58 | 9.00 to -10.00 | OEM: C, ST, CL, P |
| 9941 UL w Nosepad | 55-22-130 | 6BC 65 | 64 | SS | Y | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Airborne | 57-21-145 | 4BC 65 | 74.7 | SS | Y | Polycarbonate | Glasses | ✓ | ✓ | 62 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Centrist | 55-17-140 | 4BC 65 | 88 | Incl | Y | Nylon | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Circuit | 60-18-135 | 8BC 73 | 82.1 | SS/SSA | Y | Nylon | Glasses | ✓ | - | - | - | OEM: C, ST, CL, P |
| Mako | 65-15-130 | 8BC 73 | 69.3 | Incl | Y | Polycarbonate | Glasses | ✓ | - | - | - | OEM: C, ST, CL, P |
| Onyx | 68-19-135 | 8BC 73 | 71.7 | N/A | Y | Polycarbonate | Glasses | ✓ | - | - | - | OEM: C, ST, CL, P |
| Razer | 58-17-135 | 6BC 65 | 62.5 | Incl | Y | Polycarbonate | Glasses | ✓ | ✓ | 59 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Kids | | | | | | | | | | | | |
| Tugga Baby | 46-18 | 4BC 65 | 42 | N/A | Y | Polycarbonate | Glasses | ✓ | - | - | - | OEM: C, ST, CL, P |
| Abaco Brand | | | | | | | | | | | | |
| Dockside | 58-19-135 | 6BC 65 | 72 | N/A | Y | Polycarbonate | Glasses | ✓ | ✓ | 60 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Kai | 56-18-135 | 4BC 65 | 65 | N/A | Y | Polycarbonate | Glasses | ✓ | ✓ | 60 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Jesse | 56-18-142 | 6BC 65 | 68 | N/A | Y | Polycarbonate | Glasses | ✓ | ✓ | 60 | 8.00 to -9.00 | OEM: C, ST, CL, P |
| Waikiki | 54-21-140 | 4BC 65 | 68 | N/A | Y | Polycarbonate | Glasses | ✓ | ✓ | 60 | 8.00 to -9.00 | OEM: C, ST, CL, P |

Included Accessories:

OEM = OEM accessories are provided by Protech or customer designed accessories.

BR = Branded accessories means some or all of the included accessories are supplied by the brand.

Accessory Types:

C = Case
P = Pouch
ST = Straps
CL = Cloth

| | | | | | | | | | Prescription (Rx) Options | | | | |
|----------------------|-------------|------------|------------|--------------|----|------------------|------------|-------|---------------------------|-------------|---------------|------------------------|--|
| Model | Size | Base Curve | Weight (g) | Side Shields | CE | Material | Frame Type | Plano | Rx | Min. PRO PD | Powers | Included Accessories | |
| Costa Brand | | | | | | | | | | | | | |
| Aleta | 54-19-136 | 6BC 65 | - | N/A | N | Bioplastics | Glasses | ✓ | ✓ | - | - | - | |
| Aransas | 58-17-145 | 6BC 65 | - | N/A | N | Acetate | Glasses | ✓ | ✓ | - | - | - | |
| Blackfin | 62-17-120 | 6BC 65 | - | N/A | N | Bioplastics | Glasses | ✓ | ✓ | - | - | - | |
| Bimini Road 330 | 55-16-145 | 4BC 65 | - | SSA | N | Metal | Glasses | ✓ | ✓ | - | - | - | |
| Bimini Road 320 | 53-17-145 | 4BC 65 | - | SSA | N | Metal | Glasses | ✓ | ✓ | - | - | - | |
| Fernandina | 57-15-132 | 4BC 65 | - | SSA | N | Metal, Acetate | Glasses | ✓ | ✓ | - | - | - | |
| Ocean Ridge 400 | 57-16-140 | 4BC 65 | - | SSA | N | Bioplastics | Glasses | ✓ | ✓ | - | - | - | |
| Seadrift | 60-15-123 | 8BC 73 | - | N/A | N | Bioplastics | Glasses | ✓ | - | - | - | - | |
| Sullivan | 53-19-140 | 6BC 65 | - | SSA | N | Acetate | Glasses | ✓ | ✓ | - | - | - | |
| Tuna Alley | 61.5-11-125 | 8BC 73 | - | N/A | N | Nylon | Glasses | ✓ | - | - | - | - | |
| Tybee | 55-20-140 | 6BC 65 | - | SSA | N | Acetate | Glasses | ✓ | ✓ | - | - | - | |
| Nike Brand | | | | | | | | | | | | | |
| Brazen | 60-16-140 | 6BC 65 | 75 | SS | Y | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | BR: P ; OEM: C, ST, CL | |
| Brazen Fury | 60-17-130 | 6BC 65 | - | SS/SSA | N | Bioplastics | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | BR: P ; OEM: C, ST, CL | |
| Maverick | 59-15-145 | 6BC 65 | 70.3 | SSA | N | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | BR: P ; OEM: C, ST, CL | |
| Premier | 60-13-135 | 6BC 65 | 67.8 | N/A | N | Nylon | Glasses | ✓ | ✓ | 64 | 8.00 to -9.00 | BR: P ; OEM: C, ST, CL | |
| Rabid | 63-14-145 | 8BC 73 | 63 | SSA | Y | Nylon | Glasses | ✓ | - | - | - | BR: P ; OEM: C, ST, CL | |
| 7118 | 55-14-140 | 6BC 65 | 71.5 | SSA | N | Bioplastics | Glasses | ✓ | ✓ | 55 | 8.00 to -9.00 | BR: P ; OEM: C, ST, CL | |
| Liberty Brand | | | | | | | | | | | | | |
| Maxx 30 (small) | 53-17-130 | 6BC 65 | 60 | SS | Y | Polycarbonate | Glasses | ✓ | ✓ | 56 | - | OEM: C, ST, CL, P | |
| Maxx 30 (regular) | 55-20-130 | 6BC 65 | 69 | SS | Y | Polycarbonate | Glasses | ✓ | ✓ | 56 | - | OEM: C, ST, CL, P | |
| Maxx 31 (small) | 53-17-130 | 6BC 65 | 60 | SS | Y | Polycarbonate | Glasses | ✓ | ✓ | 56 | - | OEM: C, ST, CL, P | |
| Maxx 31 (regular) | 55-20-130 | 6BC 65 | 64 | SS | Y | Polycarbonate | Glasses | ✓ | ✓ | 56 | - | OEM: C, ST, CL, P | |
| Oakley Brand | | | | | | | | | | | | | |
| Crosslink Zero | 58-16-138 | 6BC 65 | 61.6 | N/A | Y | O Matter | Glasses | ✓ | ✓ | 56 | 8.00 to -9.00 | BR: C, ST, CL | |
| Five Squared | 54-20-135 | 8BC 73 | 71 | N/A | Y | O Matter | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| Holbrook | 56-18-135 | 6BC 65 | 66.5 | N/A | N | O Matter | Glasses | ✓ | ✓ | 60 | 8.00 to -9.00 | BR: C, ST, CL | |
| Jupiter Squared | 56-18-135 | 6BC 65 | 74.1 | N/A | Y | O Matter | Glasses | ✓ | ✓ | 66 | 8.00 to -9.00 | BR: P ; OEM: C, ST, CL | |
| Litebeam | 55-18-140 | 4BC 65 | 59.6 | N/A | N | O Matter | Glasses | ✓ | ✓ | 60 | 8.00 to -9.00 | BR: C, ST, CL | |
| Straightlink | 61-19-130 | 8BC 73 | 77.2 | N/A | Y | O Matter | Glasses | ✓ | - | - | - | OEM: C, ST, CL, P | |
| Trillbe X | 52-18-140 | 4BC 65 | 63.9 | SSA | Y | O Matter | Glasses | ✓ | ✓ | 58 | 8.00 to -9.00 | BR: C, ST, CL | |
| Ray ban Brand | | | | | | | | | | | | | |
| Daddy-O 2016 | 59-17-125 | 6BC 65 | 61 | SS | Y | Nylon | Glasses | ✓ | ✓ | 65 | 8.00 to -9.00 | BR: C ; OEM: ST, CL, P | |
| 4105 | 50-23-140 | 6BC 65 | 62 | N/A | Y | Metal | Glasses | ✓ | ✓ | 65 | 8.00 to -9.00 | BR: C, ST, CL | |
| 4165F Justin | 55-18-145 | 6BC 65 | 73 | N/A | Y | Nylon | Glasses | ✓ | ✓ | 65 | 8.00 to -9.00 | BR: C, ST, CL | |
| 4181 | 57-16-145 | 6BC 65 | 71.4 | SSA | N | Injected Plastic | Glasses | ✓ | ✓ | 65 | 8.00 to -9.00 | BR: C, ST, CL | |
| 7047 | 54-17-140 | 4BC 65 | 63 | SSA | N | Injected Plastic | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| Wiley X Brand | | | | | | | | | | | | | |
| Airrage | 61-18-124 | 8BC 73 | 67.5 | N/A | N | Polycarbonate | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| Boss | 68-18-125 | 8BC 73 | 84.5 | N/A | Y | Polycarbonate | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| Brick | 63-18-120 | 8BC 73 | 74 | N/A | Y | Polycarbonate | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| Contour | 54-17-140 | 4BC 65 | 62.3 | Incl | Y | Nylon | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| P-17 | 60-18-125 | 8BC 73 | 68 | SSA | Y | Polycarbonate | Glasses | ✓ | - | - | - | BR: C, ST, CL | |
| Face Shields / Masks | | | | | | | | | | | | | |
| 450FPM | Adjustable | - | 350g | N/A | Y | Plastic/Acrylic | Mask | ✓ | - | - | - | Custom Case | |
| 400PSM | Adjustable | - | 425g | N/A | Y | Plastic/Acrylic | Mask | ✓ | - | - | - | Custom Case | |

Included Accessories:

OEM = OEM accessories are provided by Protech or customer designed accessories.

BR = Branded accessories means some or all of the included accessories are supplied by the brand.

Accessory Types:

C = Case
P = Pouch
ST = Straps
CL = Cloth



CONTACT US

Telephone: 561.627.9769

Email (USA): sales@protechmed.com

Email (Intl.): global@protechmed.com

Email (Lat. Am.): latam@protechmed.com

www.protechmed.com

Protech Medical
1360 North Killian Dr., Unit 2
Lake Park, FL 33403