



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

NIKE 7118

MODEL: [RE-N-7118]

DESCRIPTION:

The Nike 7118 lead glasses are a retro-chic, gender-neutral frame made of lightweight, bio-plastic. Featuring rubberized nose pads and temple tips, the 7118 protective eyewear are a comfortable "all-day-wear" option.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Injected bio-plastic

Weight: 71.5g
Frame size: 55-14-140
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 55
Base Curve: 6

FEATURES:

Rubberized nose pads and temple tips for added grip and stability, lightweight frame

COLORS:

Matte Black/Blue (MBLK), Matte Midnight Navy (MNAVY)



MNAVY



MBLK

