



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Website: protechmed.com
Phone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): LatAm@protechmed.com

OAKLEY JUPITER SQUARED



MODEL: [OA-JUP]

DESCRIPTION:

A customer favorite, the sporty Jupiter Squared is attitude with refinement. The handles of switchblade knives inspired the earstems and their metal bands and rivets. The Matte Black frames of the Flag Collection are highlighted with American flag emblems.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Durable O Matter™

Weight: 74.1g
Frame size: 56-18-135
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 66
Base Curve: 8

FEATURES:

Wraparound coverage for optimal protection. Rubberized ear socks for extra comfort and grip. American flag on Oakley "O" logo.

COLORS:

Matte Black (MBLK)



MBLK

